
D 720 HELIOS

Blood Heat Exchanger

dideco

INSTRUCTIONS FOR USE

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MODE D'EMPLOI

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A. DESCRIPTION

The D 720 HELIOS is a heat exchanger for perfusion liquids. The devices are single use, non-toxic, non pyrogenic, supplied **STERILE** and packaged singly. Sterilised by ethylene oxide.

The level of ethylene oxide residuals in the device is within the limits established by national regulations in the country of use.

The device is available in the following versions:

[A] D 720A Adult

[B] D 720P Paediatric

[C] D 720C Cardioplegia

B. TECHNICAL FEATURES

Features	D 720A	D 720P	D 720C
Material	Stainless steel		
Type	Pleated steel		
Heat exchange area	0.22 sqm	0.09 sqm	0.09 sqm
Prime static volume	80 ml	43 ml	60 ml
Connections	(9.3 mm)	(6.35 mm)	(6.35 mm)
Inlet	3/8"	1/4"	1/4"
Outlet	3/8"	1/4"	

C. INTENDED USE

The D 720 HELIOS is intended for use in extracorporeal circulation for the monitoring and adjustment of perfusate.

Blood to be treated must contain anticoagulant.

The adult and paediatric versions are intended to monitor and adjust the temperature during local perfusion.

The cardioplegia version is to be used for the perfusion of blood and/or crystalloid and/or cardioplegia fluid. It is composed of a heat exchanger and an air purging section.

All versions have a probe holder for inlet and outlet temperature detection. The D 720A is an ADULT heat exchanger to be used with flows up to 8 l/min. The D 720P is a PAEDIATRIC heat exchanger to be used with flows up to 4 l/min.

The D 720C is a heat exchanger for CARDIOPLEGIA. The D 720 should not be used longer than 6 hours. Contact with blood for longer periods is not advised. The D 720 HELIOS should be used in combination with medical devices listed in section (Medical devices for use with the D 720 HELIOS).

D. SAFETY INFORMATION

Information intended to attract the attention of the user to potentially dangerous situations and to ensure correct and safe use of the device is indicated in the text in the following way:

WARNING

Indicates serious adverse reactions and potential safety hazards for practitioner and/or patient that can occur in the proper use or misuse of the device, and also the limitations of use and the measures to be adopted in such cases.

CAUTION

Indicates any special care to be exercised by a practitioner for the safe and effective use of the device.

EXPLANATION OF THE SYMBOLS USED ON THE LABELS

For single use only (Do not reuse)

Batch code (number)
(reference for product traceability)

Use by (Expiry date)

Date of manufacture

Sterile - Ethylene oxide
sterilised

Non Pyrogenic

Latex free

Warning: Do not resterilize.

Contents sterile only if package is not opened, damaged or broken

REF

Catalogue (code) number

Attention, see instruction for use

This way up

Fragile; handle with care

ea Quantity

Keep away from heat

Keep dry

The following is general safety information with the aim of advising the operator in preparing to use the device.

Also, specific safety information is given in the instructions for use at locations in the text where that information is relevant for correct operation.

WARNING

- The device must be used in accordance with the instructions for use provided in this manual.
- The device is intended to be used by professionally trained personnel.
- SORIN GROUP ITALIA is not responsible for problems arising from inexperience or improper use.
- **FRAGILE**, handle with care.
- Keep dry. Store at room temperature.
- For single use and for single-patient use only. During use the device is in contact with human blood, body fluids, liquids or gases for the purpose of eventual infusion, administration or introduction into the body, and due to its specific design it cannot be fully cleaned and disinfected after use. Therefore, reuse on other patients might cause cross-contamination, infection and sepsis. In addition, the reuse increases the probability of product failure (integrity, functionality and clinical effectiveness).

- Always give and maintain a correct dose and accurate monitoring of the anticoagulant before, and during use with blood.
- The device must not undergo any further processing.
- Do not resterilise.
- After use, dispose of the device in accordance with applicable regulations in force in the country of use.
- The device must only be used if sterile. In the event of the device being supplied NOT STERILE (in which case, this is indicated by the wording "NOT STERILE" on the pack) contact DIDE-CO or an authorised dealer to agree on a sterilisation method.
- For further information and/or in case of complaint contact SORIN GROUP ITALIA or the authorised local representative.

CAUTION

Federal law (U.S.A.) restricts this device to sale by or on the order of a physician.

E. SET-UP

1) POSITION THE HOLDER

Position the HELIOS D 611 holder on the pump structure by means of the clamp.

2) FIX TO THE HOLDER

WARNING

- Sterility is guaranteed only if the sterile packaging is not wet, opened, damaged or broken. Do not use the device if sterility cannot be guaranteed.
- Check the expiry date on the label attached. Do not use the device after the date shown.
- The device must be used immediately after opening the sterile packaging.
- The device must be handled aseptically.

Remove the device from the sterile packaging.

WARNING

- Carry out a visual inspection and carefully check the device before use. Transport and/or storage conditions other than those prescribed may have caused damage to the device.
- Do not use solvents such as alcohol, ether, acetone, etc.: as contact may cause damage to the device.
- Do not allow halogenated liquids such as

Halothane and Fluothane to come into contact with the polycarbonate housing of the device. This could cause damage which may compromise the integrity and proper function of the device.

Remove the protective caps from the heat exchanger water connectors.

Insert the HELIOS in the holder support clamp.

Close the clamp by means of the side knob.

3) THERMOCIRCULATOR SET UP

Connect the water tubes to the holder by means of the female Hansen connectors SORIN GROUP ITALIA code 9028 (D 720A/P - ref. 2, D 720C - ref.3).

WARNING

- The use of different connectors from those indicated may cause resistance inside the circuit and reduce the efficiency of the heat exchanger.
- Do not occlude the hole on the heat exchanger bottom cover as it is the outlet of the safety chamber which prevents contamination in the event of a seal failure.
- Water pressure in the device must not exceed 300 kPa (3 bar / 44 psi).
- Pressure in the blood compartment of the heat exchanger must not exceed 100 KPa (1 bar / 14 psi).
- The water temperature at the heat exchanger inlet must not exceed 42°C (108°F).

4) CHECK THE HEAT EXCHANGER

Check the heat exchanger by recirculating water inside the heat exchanger for a few minutes.

The integrity of the housing is guaranteed if there are absolutely no leaks from the water compartment or from the hole of the safety channel.

5) HEAT EXCHANGER CONNECTION

WARNING

All connections downstream of the pump must be secured by means of ties.

INLET LINE: connect a 3/8" (1/4") line to the inlet connector (ref.1) of the D 720A (D 720P/C).

OUTLET LINE: connect a 3/8" (1/4") line to the outlet connector (ref.2) of the D 720A (D 720P/C).

Cardioplegia version

PURGE LINE: remove the protective cap of luer connector (ref. 6) on the bubble trap head and connect the male connector of a purge line.

6) CONNECT THE TEMPERATURE PROBES

Inlet/outlet blood temperature probes site are beside the inlet/outlet connectors. SORIN GROUP ITALIA temperature probes code no. is 9026.

F. PRIMING PROCEDURE

1) CONNECTION TO THE EXTRACORPOREAL CIRCUIT

Adult and Paediatric versions

Connect the D 720 A/P to the extracorporeal circuit.

Cardioplegia version

Connect D720 C inlet to the coronary outlet (if available) or to an oxygenator arterial outlet by means of a dedicated pump. Connect D720C outlet to the cardioplegia cannulae.

2) CIRCUIT PRIMING

Prime the circuit.

All versions

Turn the heat exchanger upside down and start priming at a 20 - 30 ml/min flow. As soon as the liquid flows out of the outlet connector clamp the outlet tube and turn the heat exchanger to its original position.

Cardioplegia version

Open the purge line place on the bubble trap head.

3) PURGE THE AIR

Tap the heat exchanger and open the outlet line .

Adult and paediatric versions

Direct air towards the outlet to facilitate the purging and keep the outlet upwards. This procedure can last for some minutes.

The procedure is successful when air bubbles stop exiting the blood outlet connector

Cardioplegia version

During operation purge the air from the bubble trap by opening the purge line from time to time.

CAUTION

The most effective bubble purging technique demands for a high flow recirculation of and therefore a circuit which allows for it. Keep the heat exchanger upside down during recirculation. It is necessary to tap the device until bubbles are completely evacuated.

WARNING

The ACT (Activated Coagulation Time) must always be longer than or equal to 480 seconds

in order to ensure adequate anticoagulation of extracorporeal circuit.

G. HEAT EXCHANGER REPLACEMENT

A spare heat exchanger must always be available during perfusion. After 6 hours of use with blood or if particular situations occur, which may lead the person responsible for perfusion to determine that the safety of the patient may be compromised (insufficient performance, leaks from the lower cover hole of D720, etc.), proceed as follows for the replacement of the device.

CAUTION

Use a sterile procedure during the entire replacement phase.

- 1) Stop the pump
- 2) Place two clamps on the inlet and outlet lines (at a 5 cm distance).
- 3) Switch the thermocirculator off, clamp water lines and remove them.
- 4) Remove the purge line (only D 720 C) and all monitoring lines.
- 5) Cut inlet and outlet lines in a point between the two clamps.
- 6) Remove the HELIOS to be replaced from the holder.
- 7) Fix the new HELIOS to the holder. Connect all lines (i.e. inlet line, outlet line, purge line for the D720 C).
- 8) Remove the clamps from the water line and switch the thermocirculator on; check the new HELIOS for water leaks.
- 9) Prime the new HELIOS and evacuate any air bubble following the priming procedure instructions.
- 10) Check the connections and seal them by means of ties.
- 11) Restore the perfusion.

H. MEDICAL DEVICES FOR USE WITH THE D 720 HELIOS

Cardioplegia version

The device must be used in combination with a purge line with a luer lock end to be connected to a luer connector on the bubble trap head. I

All versions

All tubing used to make the circuit connections must be of a diameter which is compatible with the dimensions of the connectors on the device (3/8"

or 1/4"). Temperature controls must be carried out using SORIN GROUP ITALIA probes, code 9026 or compatible to YSI Series 400. Any heating/cooling system (thermocirculator) may be used, provided that the connectors to the water distributor holder are of the Hansen type (SORIN GROUP ITALIA code 9028). Currently SORIN GROUP ITALIA is not aware of any contraindications to the use of the device with occlusive or non-occlusive peristaltic pumps or with centrifugal pumps. The use of other types of pump must be agreed with SORIN GROUP ITALIA.

I. RETURN OF USED PRODUCTS

Should the user be dissatisfied with anything related to the quality of the product, the product distributor or the authorized local SORIN GROUP ITALIA representative should be notified.

All parameters considered critical by the user must be reported with particular care and urgency. The following is the minimum information that should be provided:

- Detailed description of the event and, if pertinent, the conditions of the patient;
- Identification of the product involved;
- Lot number of the product involved;
- Availability of the product involved;
- All the indications the user considers useful in order to understand the origin of the elements of dissatisfaction.

SORIN GROUP ITALIA reserves the right to authorize, if necessary, recall of the product involved in the notification for assessment. If the product to be returned is contaminated, it must be treated, packed and handled in conformity with the provisions of the legislation in force in the country where the product was used.

CAUTION

It is the responsibility of the health care institution to adequately prepare and identify the product for return shipment. Do not return products that have been exposed to blood borne infectious diseases.

ONLY for US customers

If for any reason the product must be returned to the manufacturer, a returned good authorisation (RGA) number is required from Sorin Group USA cardiovascular Inc. prior to shipping.

If the product has been in contact with blood or blood fluids, it must be thoroughly cleaned and disinfected before packing. It should be shipped in either the original carton or an equivalent carton to prevent damage during shipment, and it should be properly labelled with an RGA number and an indication of the biohazardous nature of the content in the shipment.

The shipping address for returned goods in the US is:

Sorin Group USA, Inc.
Returned CV Products
14401 West 65th Way
Arvada, CO 80004-3599
FAX (800) 323 4031.

L. LIMITED WARRANTY

This Limited Warranty is in addition to any statutory rights of the Purchaser pursuant to applicable law.

SORIN GROUP ITALIA warrants that all reasonable care has been taken in the manufacture of this medical device, as required by the nature of the device and the use for which the device is intended. SORIN GROUP ITALIA warrants that the medical device is capable of functioning as indicated in the current instructions for use when used in accordance with them by a qualified user and before any expiry date indicated on the packaging. However, SORIN GROUP ITALIA cannot guarantee that the user will use the device correctly, nor that the incorrect diagnosis or therapy and/or that the particular physical and biological characteristics of an individual patient, do not affect the performance and effectiveness of the device with damaging consequences for the patient, even though the specified instructions for use have been respected.

SORIN GROUP ITALIA, whilst emphasizing the need to adhere strictly to the instructions for use and to adopt all the precautions necessary for the correct use of the device, cannot assume any responsibility for any loss, damage, expense, incidents or consequences arising directly or indirectly from the improper use of this device.

SORIN GROUP ITALIA undertakes to replace the medical device in the event that it is defective at the time of placing on the market or whilst being shipped by SORIN GROUP ITALIA up to the time of delivery to the final user unless such defect has been

caused by mishandling by the purchaser. The above replaces all other warranties explicit or implicit, written or verbal, including warranties of merchantability and fitness for purpose. No person, including any representative, agent, dealer, distributor or intermediary of SORIN GROUP ITALIA or any other industrial or commercial organization is authorized to make any representation or warranty concerning this medical device except as expressly stated herein.

SORIN GROUP ITALIA disclaims any warranty of merchantability and any warranty of fitness for purpose with regard to this product other than what is expressly stated herein.

The purchaser undertakes to comply with the terms of this Limited Warranty and in particular agrees, in the event of a dispute or litigation with SORIN GROUP ITALIA, not to make claims based on alleged or proven changes or alterations made to this Limited Warranty by any representative, agent, dealer, distributor or other intermediary.

The existing relations between the parties to the contract (also in the case that it is not drawn up in writing) to whom this Warranty is given as well as every dispute related to it or in any way connected to it as well as anything related to it or any dispute concerning this Warranty, its interpretation and execution, nothing excluded and/or reserved, are regulated exclusively by the Italian law and jurisdiction. The court chosen is the Court of Modena (Italy).

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

D 720 HELIOS

Adult and Pediatric size Blood Heat Exchanger.




- 1 - Blood inlet and outlet connectors
3/8" I.D. (Adult) and 1/4" I.D. (Pediatric)
- 2 - Water connectors "Hansen" type.
- 3 - Blood temperature probe sites.

D 720 HELIOS

Cardioplegia Heat Exchanger.

- 1 - Blood inlet connector 1/4" I.D.
- 2 - Blood outlet connector 1/4" I.D.
- 3 - Water connectors "Hansen" type.
- 4 - Temperature probe site: blood in.
- 5 - Temperature probe site: blood out.
- 6 - Purging female luer lock.

This medical device bears the **CE** marking according to the European Council Directive MDD 93/42/EEC.
Further information is available from Manufacturer (contact Sorin Group Italia's local Representative or directly Sorin Group Italia's RA & QA department).

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