BRAT2[®] Autologous Blood Recovery Suction Assemblies

Catalog Numbers 007005001 007005400

INSTRUCTIONS FOR USE

DESCRIPTION

The BRAT[®] Suction Assemblies consist of a suction tip connector bonded to a dual lumen tube assembly. The small lumen of the dual lumen tube supplies anticoagulant solution to the salvaged blood. Flow control for the anticoagulant solution is achieved with the use of a roller clamp and universal drip chamber spike. The large lumen tube connects to the 1/4" inlet port of a blood collection reservoir and facilitates the transfer of salvaged blood mixed with anticoagulant from the sucker tip to the blood collection reservoir.

Contents of this set have been sterilized by ethylene oxide gas and have nonpyrogenic fluid pathways.

INDICATIONS FOR USE

BRAT Suction Assemblies are intended for use in recovery and/or processing of autologous blood.

CONTRAINDICATIONS

There are no known contraindications for this product. However, the use of blood from this device for reinfusion may be contraindicated, e.g., in the presence of certain potentially harmful substances (refer to the BRAT2[®] Operator's Manual for more detail). The responsibility for the use of the device in all cases belongs solely to the physician in charge.

WARNINGS

Carefully observe this set for leaks before and during use. Leakage may result in loss of sterility or loss of blood and/or fluid. If leakage is observed before or during use, replace leaking component or retighten leaking connection as appropriate.

PRECAUTIONS

1. Carefully read these Instructions for Use before using this product. Refer to the BRAT2 Operator's Manual for complete instructions.
2. Due to the possibility of operator exposure to bloodborne pathogens (such as HIV, hepatitis viruses, bacteria, Cytomegalovirus, etc.) when handling extracorporeal blood circuits, adequate precautions should be taken at all times to prevent the exposure to and transmission of such agents.
3. Discard after single use. Do not resterilize any part of this set, including packaging and wrapping material.
4. Do not use this set if it is damaged, if the sterile package is damaged or opened, or if the protective

caps are not in place. Check for loose connections on preassembled kits.

5. This product is intended for use by trained personnel only.
6. Use proper aseptic technique while handling this set.
7. If using a spike with capped vent, observe the following: When using a rigid IV bottle, the cap on the spike vent should be open to allow air to displace fluid as it drains from the bottle.
8. Federal law (U.S.A.) restricts this device to sale by or on the order of a physician.

PROCEDURES FOR USE

007005001

1. Using aseptic technique, prepare a bag or bottle of anticoagulant solution by mixing 30,000 units of heparin with a liter of sterile, isotonic normal saline (injection) solution or lesser amounts in same proportion. Premixed ACD-A or CPD solutions in bags or bottles can be used instead of the heparin solution.
2. Using aseptic technique, remove BRAT Suction Assembly from its package, remove outer wrap, and pass the inner sterile wrapped assembly to the sterile field.
3. At the sterile field, unwrap assembly, remove red protective cover on BRAT Suction Assembly, and attach connector to suction tip. Pass other end of BRAT Suction Assembly from the sterile field.
4. Remove blue protective cover on BRAT Suction Assembly and attach to a 1/4" defoamed inlet port of a blood collection reservoir.
5. Clamp the solution line. Using aseptic technique, spike the container of anticoagulant solution and hang it a minimum of two feet above the surgical field. Fill the drip chamber 1/4 full. Prior to use, prime the tubing and reservoir with sufficient solution to wet the blood contact surfaces (100-200ml). Assure suction is on during the priming operation.
6. Anticoagulation of salvaged blood should be carried out by the standard protocol within the institution using the Suction Assembly. In the absence of such a protocol, the following guidelines are suggested:

Heparin. Adjust and monitor the flow rate during use such that one part heparin solution is mixed with four parts shed blood. A rate of one drop per second (20 drops/ml) is sufficient to anticoagulate an end blood product in the reservoir of approximately one liter per hour.

Citrate. Adjust and monitor the flow rate during use such that one part of ACD-A or CPD solution is mixed

with 5 to 10 parts blood. A rate of 2 drops every 3 seconds (20 drops/ml) is sufficient to anticoagulate an end blood product in the reservoir of approximately one liter per hour.

007005400

- Using aseptic technique, prepare a bag or bottle of anticoagulant solution by mixing 30,000 units of heparin with a liter of sterile, isotonic normal saline (injection) solution or lesser amounts in same proportion. Premixed ACD-A or CPD solutions in bags or bottles can be used instead of the heparin solution.
- Using aseptic technique, remove BRAT Suction Assembly from its package, remove outer wrap, and pass the inner sterile wrapped assembly to the sterile field.
- At the sterile field, unwrap assembly, remove red protective cover on BRAT Suction Assembly, and attach connector to suction tip. Pass other end of BRAT Suction Assembly from the sterile field.
- Remove blue protective cover on BRAT Suction Assembly and attach to a 1/4" defoamed inlet port of a Sorin Group Italia blood collection reservoir. Turn the Quick-Lock Connector to secure the line to the inlet.
- Clamp the solution line. Using aseptic technique, spike the container of anticoagulant solution and hang it a minimum of two feet above the surgical field. Fill the drip chamber 1/4 full. Prior to use, prime the tubing and reservoir with sufficient solution to wet the blood contact surfaces (100-200ml). Assure suction is on during the priming operation.
- Anticoagulation of salvaged blood should be carried out by the standard protocol within the institution using the Suction Assembly. In the absence of such a protocol, the following guidelines are suggested:

Heparin. Adjust and monitor the flow rate during use such that one part heparin solution is mixed with four parts shed blood. A rate of one drop per second (20 drops/ml) is sufficient to anticoagulate an end blood product in the reservoir of approximately one liter per hour.

Citrate. Adjust and monitor the flow rate during use such that one part of ACD-A or CPD solution is mixed with 5 to 10 parts blood. A rate of 2 drops every 3 seconds (20 drops/ml) is sufficient to anticoagulate an end blood product in the reservoir of approximately one liter per hour.

REFERENCE ARTICLES

- "Minimizing Blood Transfusions During Abdominal Aortic Surgery: Recent Advances in Rapid Autotransfusion," Hallett J.W., et al, *Journal of Vascular Surgery*, 5(4): 601, 1987.
- "Panel IV: Dosage Regimens for Citrate Anticoagulants," Hester, J.P., et al, *Journal of Clinical Apheresis*, 1: 149-157, 1983.
- "Evaluation of a New Blood Autotransfusion Device," Toomasian, J.M. et al, *Journal of Thoracic and Cardiovascular Surgery*, 92: 936-943, 1986.

SPECIFICATIONS

Dual Lumen Tubing		
Large	9.5mm O.D.	6.4mm I.D.
Small	4.3mm O.D.	3.0mm I.D.
Roller Clamp	Control flow rate of 4.3mm O.D. tube	
Aspiration Connector	Open end 6.6mm	

RETURN OF USED PRODUCT

Should the user be dissatisfied with anything related to the quality of the product, the product distributor or the authorized local SORIN GROUP ITALIA representative should be notified.

All parameters considered critical by the user must be reported with particular care and urgency. The following is the minimum information that should be provided:

- Detailed description of the event and, if pertinent, the conditions of the patient;
- Identification of the product involved;
- Lot number of the product involved;
- Availability of the product involved;
- All the indications the user considers useful in order to understand the origin of the elements of dissatisfaction.

SORIN GROUP ITALIA reserves the right to authorize, if necessary, recall of the product involved in the notification for assessment. If the product to be returned is contaminated, it must be treated, packed and handled in conformity with the provisions of the legislation in force in the country where the product was used.

CAUTION

It is the responsibility of the health care institution to adequately prepare and identify the product for return shipment. Do not return products that have been exposed to blood borne infectious diseases.

ONLY for US customers

If for any reason the product must be returned to the manufacturer, a returned good authorisation (RGA) number is required from Sorin Group USA, Inc. prior to shipping. If the product has been in contact with blood or blood fluids, it must be thoroughly cleaned and disinfected before packing. It should be shipped in either the original carton or an equivalent carton to prevent damage during shipment, and it should be properly labelled with an RGA number and an indication of the biohazardous nature of the content in the shipment.

CAUTION

It is the responsibility of the health care institution to adequately prepare and identify the product for return shipment. Do not return products that have been exposed to blood borne infectious diseases.

The shipping address for returned goods in the US is:

SORIN GROUP USA, Inc.
 Returned CV Products
 14401 West 65th Way
 Arvada, CO 80004-3599
 FAX (800) 323 4031

LIMITED WARRANTY

This Limited Warranty is in addition to any statutory rights of the Purchaser pursuant to applicable law.

SORIN GROUP ITALIA warrants that all reasonable care has been taken in the manufacture of this medical device, as required by the nature of the device and the use for which the device is intended.

SORIN GROUP ITALIA warrants that the medical device is capable of functioning as indicated in the current instructions for use when used in accordance with them by a qualified user and before any expiry date indicated on the packaging. However, SORIN GROUP ITALIA cannot guarantee that the user will use the device correctly, nor that the incorrect diagnosis or therapy and/or that the particular physical and biological haematological characteristics of an individual patient, do not affect the performance and effectiveness of the device with damaging consequences for the patient, even though the specified instructions for use have been respected.

SORIN GROUP ITALIA, whilst emphasizing the need to adhere strictly to the instructions for use and to adopt all the precautions necessary for the correct use of the device, cannot assume any responsibility for any loss, damage, expense, incidents or consequences arising directly or indirectly from the improper use of this device.

SORIN GROUP ITALIA undertakes to replace the medical device in the event that it is defective at the time of placing on the market or whilst being shipped by SORIN GROUP ITALIA up to the time of delivery to the final user unless such defect has been caused by mishandling by the purchaser.

The above replaces all other warranties explicit or implicit, written or verbal, including warranties of merchantability and fitness for purpose. No person, including any representative, agent, dealer, distributor or intermediary of SORIN GROUP ITALIA or any other industrial or commercial organization is authorized to make any representation or warranty concerning this medical device except as expressly stated herein. SORIN GROUP ITALIA disclaims any warranty of merchantability and any warranty of fitness for purpose with regard to this product other than what is expressly stated herein. The purchaser undertakes to comply with the terms of this Limited Warranty and in particular agrees, in the event of a dispute or litigation with SORIN GROUP ITALIA, not to make claims based on alleged or proven changes or alterations made to this Limited Warranty by any representative, agent, dealer, distributor or other intermediary.

The existing relations between the parties to the contract (also in the case that it is not drawn up in writing) to whom this Warranty is given as well as every dispute related to it or in any way connected to it as well as anything related to it or any dispute concerning this Warranty, its interpretation and execution, nothing excluded and/or reserved, are regulated exclusively by the Italian law and jurisdiction. The court chosen is the Court of Modena (Italy).

DEFINITION OF SYMBOLS

(as used in product labeling)

	FOR SINGLE USE ONLY (DO NOT REUSE)
	BATCH CODE (NUMBER) (REFERENCE FOR PRODUCT TRACEABILITY)
	USE BY (EXPIRY DATE)
	DATE OF MANUFACTURE
	STERILE - ETHYLENE OXIDE STERILISED
	NON PYROGENIC
	CONTAINS PHTHALATE
	LATEX FREE
Rx ONLY	CAUTION: FEDERAL LAW (USA) RESTRICT THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN.
REF	CATALOGUE (CODE) NUMBER
	ATTENTION, SEE INSTRUCTION FOR USE
	THIS WAY UP
	FRAGILE; HANDLE WITH CARE
	KEEP AWAY FROM HEAT
EA	QUANTITY
	KEEP DRY
	CONTENTS STERILE ONLY IF PACKAGE IS NOT OPENED, DAMAGED OR BROKEN

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LITERATURHINWEISE

"Minimizing Blood Transfusions During Abdominal Aortic Surgery: Recent Advances in Rapid Autotransfusion," Hallett J.W., et al, *Journal of Vascular Surgery*, 5(4): 601, 1987.

"Panel IV: Dosage Regimens for Citrate Anticoagulants," Hester, J.P., et al, *Journal of Clinical Apheresis*, 1: 149-157, 1983.

"Evaluation of a New Blood Autotransfusion Device," Toomasian, J.M. et al, *Journal of Thoracic and Cardiovascular Surgery*, 92: 936-943, 1986.

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"Minimizing Blood Transfusions During Abdominal Aortic Surgery: Recent Advances in Rapid Autotransfusion," Hallett J.W., et al., *Journal of Vascular Surgery*, 5(4): 601, 1987.

"Panel IV: Dosage Regimens for Citrate Anticoagulants," Hester, J.P., et al., *Journal of Clinical Apheresis*, 1: 149-157, 1983.

"Evaluation of a New Blood Autotransfusion Device," Toomasian, J.M. et al., *Journal of Thoracic and Cardiovascular Surgery*, 92: 936-943, 1986.

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