

# Perfusion Tubing Systems

011625/06 – 02/2011

GB – ENGLISH - INSTRUCTIONS FOR USE	3	CZ – ČESKÝ JAZYK - NÁVOD K POUŽITÍ	39
IT – ITALIANO - ISTRUZIONI PER L'USO	6	PL – POLSKI - INSTRUKCJA OBSŁUGI	42
FR – FRANÇAIS - MODE D'EMPLOI	9	SK – SLOVENSKÝ JAZYK – POKYNY PRE POUŽITIE	45
DE – DEUTSCH - GEBRAUCHSANWEISUNG	12	SL - SLOVENŠČINA - NAVODILO ZA UPORABO	48
ES – ESPAÑOL - INSTRUCCIONES DE USO	15	HU – MAGYAR - HASZNÁLATI UTASÍTÁS	51
PT – PORTUGUÊS - INSTRUÇÕES PARA UTILIZAÇÃO	18	LT – LIETUVIŠKAI - NAUDOJIMO INSTRUKCIJOS	54
GR – ΕΛΛΗΝΙΚΑ - ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ	21	EE – EESTI - KASUTUSJUHEND	57
NL – NEDERLANDS - GEBRUIKSAANWIJZINGEN	24	LV – LATVISKI – LIETOŠANAS INSTRUKCIJAS	60
SE – SVENSKA - BRUKSANVISNING	27	RO – ROMÂNĂ – INSTRUCȚIUNI DE UTILIZARE	63
DK – DANSK - BRUGSANVISNING	30	BG – БЪЛГАРСКИ – ИНСТРУКЦИИ ЗА УПОТРЕБА	66
FI – SUOMALAINEN - KÄYTTÖOHJEET	33	TR - TÜRKÇE - KULLANMA TALİMATI	69
NO - NORSK - BRUKSINSTRUKSER	36		

# GB – ENGLISH - INSTRUCTIONS FOR USE

## I. CONTENTS

- I. Contents
- A. Description
- B. Intended use
- C. Safety information
- D. Set up
- E. Device replacement
- F. Return of used products
- G. Limited Warranty

## A. DESCRIPTION

This device is manufactured to user design and specifications and may consist of any or all of the following main components: oxygenators, cardiotomy reservoirs, venous reservoirs, tubing, connectors, filters, isolators, drip chambers, clamps, stopcocks, heat exchanger and bubble traps amongst others according to requirements.

Blood contact surfaces of the PTS Ph.I.S.I.O. coated have been modified to improve blood compatibility, resulting in reduced platelet adhesion on the coated surfaces.

The device is single use, non-toxic, non pyrogenic, supplied **STERILE**. Sterilised by ethylene oxide. The level of ethylene oxide residuals in the device is within the limits established by national regulations in the country of use.

## B. INTENDED USE

This tubing pack is intended for use during surgery requiring cardiopulmonary or other surgical techniques. It is intended for use together with such equipment as blood pumps, oxygenators, reservoirs, filters, heat exchangers and cannulae. The specific indications for use of individual item should be consulted.

Blood in contact with the device should contain anticoagulant.

The device should not be used longer than 6 hours. Contact with blood for longer periods is not advised.

## C. SAFETY INFORMATION

Information intended to attract the attention of the user to potentially dangerous situations and to ensure correct and safe use of the device is indicated in the text in the following way:

### ⚠ WARNING

- Indicates serious adverse reactions and potential safety hazards for the user and/or the patient that can occur in the proper use or misuse of the device, and also the limitations of use and the measures to be adopted in such cases.

### ⚠ CAUTION

- Indicates any special care to be exercised by the user for the safe and effective use of the device.

### EXPLANATION OF THE SYMBOLS USED ON THE LABELS

 For single use only (Do not reuse)

 Batch code (number)  
(reference for product traceability)

 Use by (Expiry date)

 Date of manufacture

 Manufactured by:

 Sterile - Ethylene oxide sterilised

 Non Pyrogenic

 Contains PHTHALATE

 Latex free

 Warning: Do not resterilize.

 Contents sterile only if package is not opened, damaged or broken

 Catalogue (code) number

 Attention, see instruction for use

 Quantity

 This way up

 Fragile; handle with care

 Keep away from heat

 Keep dry

The following is general safety information with the aim of advising the user when preparing to use the device.

Also, specific safety information is given in the instructions for use at locations in the text where that information is relevant for correct operation.

### ⚠ WARNING

- Only use the unit if it is **STERILE**.
- The device must be used in accordance with the instructions for use provided in this manual.
- The device is intended to be used by professionally trained personnel.
- Sorin Group Italia is not responsible for problems arising from inexperience or improper use.
- **FRAGILE**, handle with care.
- Keep dry. Store at room temperature.
- For single use and for single-patient use only. During use the device is in contact with human blood, body fluids, liquids or gases for the purpose of eventual infusion, administration or introduction into the body, and due to its specific design it cannot be fully cleaned and disinfected after use. Therefore, reuse on other patients might cause cross-contamination, infection and sepsis. In addition, the reuse increases the probability of product failure (integrity, functionality and clinical effectiveness).
- The device contains phthalates. Considering the nature of contact with the body, the limited duration of contact and the number of treatments per patient, the amount of phthalates which might be released from the device does not raise specific concerns about residual risks. Further information is available on request from Sorin Group Italia
- Always administer and maintain a correct dose and accurate monitoring of the anticoagulant before, during and after the bypass.
- The device must not undergo any further processing.
- Do not resterilise.

- After use, dispose of the device in accordance with applicable regulations in force in the country of use.
- Sterility is guaranteed only if the sterile packaging is not wet, opened, damaged or broken. Do not use the device if sterility cannot be guaranteed.
- Check the expiry date on the label attached. Do not use the device after the date shown.
- The device must be used immediately after opening the sterile packaging.
- The device must be handled aseptically.
- Carry out a visual inspection and carefully check the device before use. Transport and/or storage conditions other than those prescribed may have caused damage to the device.
- Do not use solvents such as alcohol, ether, acetone, etc.: as contact may cause damage to the device.
- Do not allow halogenated liquids such as Halothane and Fluothane to come into contact with the polycarbonate housing of the device. This could cause damage which may compromise the integrity and proper functioning of the device.
- The instructions for use of individual components should be read and understood also.
- Carefully read and follow the instructions for use of the manufacturer of the Heart Lung Machine (HLM) in use regarding tubing specifications and installation.
- Handle all tubing carefully to avoid damaging it.
- The user is responsible for any changes made either to the set or to the procedure that may affect the function of the set.
- Carefully check for leaks before and during use. Leaks may cause loss of sterility, loss of perfusate, or air emboli. Should leakage be observed before or during use, re-tighten or replace the leaking component following good perfusion practice. The use of an arterial filter or bubble trap is recommended to reduce the possibility of transmitting gaseous emboli to the patient when used in accordance with manufacturers instructions.
- If the set contains a pre-bypass filter, it must not be exposed to blood or other cellular fluids and must be removed before initiating bypass.
- The effectiveness of myocardial cooling may only be ascertained by the measurement of myocardial temperature. It is affected by variables such as, but not limited to: use and degree of whole body hypothermia, cannulation and topical cooling, surgical technique and cardiac anatomy and pathology.
- For further information and/or in case of complaint contact **SORIN GROUP ITALIA** or the authorised local representative.

#### ⚠ CAUTION

- Federal law (U.S.A.) restricts this device to sale by or on the order of a physician.

## D. SET UP

#### ⚠ WARNING

When setting-up the tubing pack take care to avoid kinking, cutting, piercing, puncturing of the tubing and other damage to any individual component part of the device.

#### ⚠ CAUTION

- Tubing used as pump segment in roller pump applications may get damaged or may lead to patient injury, if use is prolonged or installation is improper.

The following factors have been found to influence tubing life in roller pump applications:

- Type of fluid being pumped/ transported.
- Occlusion degree of tubing walls.
- Pump speed (roller impacts/minute).
- Downstream backpressure.
- Temperature of transported fluid and operating theatre.
- Wear due to friction of rollers and guides (also caused by equipment improper maintenance or wearing).
- Excessive length of tubing installed into pump raceway.
- Wall thickness and wall thickness size tolerance.

#### ⚠ WARNING

- **If this set is to be used with a roller pump, then carefully set up and monitor each of the foregoing variables during use of the tubing and observe the following:**
  - a) Inspect tubing for any damage, if none install in pump. Also monitor during use for any signs of excessive wear, cracks, leaks or air entrainment and take appropriate action.
  - b) It is recommended that pump occlusion is adjusted as per the HLM manufacturer's instructions prior to each procedure. Improper occlusion settings may cause pump segment tube failure, excessive spallation, haemolysis and/or inaccurate blood flow.
  - c) When placing tubing in a rollerpump, maintain the natural curvature of the tube and avoid twisting of the tube.
  - d) The roller pump should be maintained in accordance with the HLM manufacturer's recommendations in order to avoid increased tubing wear and possible premature failure.
  - e) The correct size or adjustment of the tube clamps should be used to avoid damage to the tube and to hold it securely.
- **If the set contains a cooling coil use only salt-free water and ice. Do not add alcohol. Reducing the ice bath temperature below 0°C may cause cold-induced injury to the patient.**
  1. Remove the set from the package using aseptic technique.
  2. When making any connections use aseptic technique.
  3. Connect luer locks or fittings finger tight. Do not over tighten, but make sure connection is secure.
  4. Position stopcock handles as appropriate.
  5. Connect the set as appropriate to oxygenator, heat exchangers, filters and other components in accordance with their instructions for use.
  6. When making any connections, ensure component parts are properly joined.
  7. Ensure any one-way valves in the set are pointing in the correct direction.
  8. Prime the circuit in accordance with the instructions relating to the oxygenator and any other components of the circuit.
  9. Inspect to make sure that all parts to be primed are free of air bubbles.
  10. Inspect to make sure that there are no leaks.
  11. Initiate cardiopulmonary bypass in accordance with oxygenator instructions and good perfusion practice.
- **Select appropriate occlusion and RPM settings that will achieve desired flow rate with the minimum stress to the tubing.**
- **Carefully re-set occlusion whenever the tubing is either moved or changed in the pump raceway.**

## E DEVICE REPLACEMENT

A spare device must always be available during perfusion. After 6 hours of use with blood or if situations occur, which lead the person responsible for perfusion to determine that the safety of the patient may be compromised, consider device replacement.

#### ⚠ CAUTION

- **Use sterile methods during all replacement procedure.**

## F. RETURN OF USED PRODUCTS

Should the user be dissatisfied with anything related to the quality of the product, the product distributor or the authorized local SORIN GROUP ITALIA representative should be notified.

All parameters considered critical by the user must be reported with particular care and urgency. The following is the minimum information that should be provided:

- Detailed description of the event and, if pertinent, the conditions of the patient;
- Identification of the product involved;
- Lot number of the product involved;
- Availability of the product involved;
- All the indications the user considers useful in order to understand the origin of the elements of dissatisfaction.

SORIN GROUP ITALIA reserves the right to authorize, if necessary, recall of the product involved in the notification for assessment. If the product to be returned is contaminated, it must be treated, packed and handled in conformity with the provisions of the legislation in force in the country where the product was used.

### ⚠ CAUTION

- **It is the responsibility of the health care institution to adequately prepare and identify the product for return shipment. Do not return products that have been exposed to blood borne infectious diseases.**

#### ONLY for US customers

If for any reason the product must be returned to the manufacturer, a returned good authorisation (RGA) number is required from Sorin Group USA, Inc. prior to shipping.

If the product has been in contact with blood or blood fluids, it must be thoroughly cleaned and disinfected before packing. It should be shipped in either the original carton or an equivalent carton to prevent damage during shipment, and it should be properly labelled with an RGA number and an indication of the biohazardous nature of the content in the shipment.

### ⚠ CAUTION

- **It is the responsibility of the health care institution to adequately prepare and identify the product for return shipment. Do not return products that have been exposed to blood borne infectious diseases.**

The shipping address for returned goods in the US is:

Sorin Group USA, Inc.  
Returned CV Products  
14401 West 65th Way  
Arvada, CO 80004-3599  
FAX (800) 323 4031

## G LIMITED WARRANTY

*This Limited Warranty is in addition to any statutory rights of the Purchaser pursuant to applicable law.*

SORIN GROUP ITALIA warrants that all reasonable care has been taken in the manufacture of this medical device, as required by the nature of the device and the use for which the device is intended.

SORIN GROUP ITALIA warrants that the medical device is capable of functioning as indicated in the current instructions for use when used in accordance with them by a qualified user and before any expiry date indicated on the packaging.

However, SORIN GROUP ITALIA cannot guarantee that the user will use the device correctly, nor that the incorrect diagnosis or therapy and/or that the particular physical and biological characteristics of an individual patient, do not affect the performance and effectiveness of the device with damaging consequences for the patient, even though the specified instructions for use have been respected.

SORIN GROUP ITALIA, whilst emphasizing the need to adhere strictly to the instructions for use and to adopt all the precautions necessary for the correct use of the device, cannot assume any responsibility for any loss, damage, expense, incidents or consequences arising directly or indirectly from the improper use of this device.

SORIN GROUP ITALIA undertakes to replace the medical device in the event that it is defective at the time of placing on the market or whilst being shipped by SORIN GROUP ITALIA up to the time of delivery to the final user unless such defect has been caused by mishandling by the purchaser.

The above replaces all other warranties explicit or implicit, written or verbal, including warranties of merchantability and fitness for purpose. No person, including any representative, agent, dealer, distributor or intermediary of SORIN GROUP ITALIA or any other industrial or commercial organization is authorized to make any representation or warranty concerning this medical device except as expressly stated herein. SORIN GROUP ITALIA disclaims any warranty of merchantability and any warranty of fitness for purpose with regard to this product other than what is expressly stated herein. The purchaser undertakes to comply with the terms of this Limited Warranty and in particular agrees, in the event of a dispute or litigation with SORIN GROUP ITALIA, not to make claims based on alleged or proven changes or alterations made to this Limited Warranty by any representative, agent, dealer, distributor or other intermediary.

The existing relations between the parties to the contract (also in the case that it is not drawn up in writing) to whom this Warranty is given as well as every dispute related to it or in any way connected to it as well as anything related to it or any dispute concerning this Warranty, its interpretation and execution, nothing excluded and/or reserved, are regulated exclusively by the Italian law and jurisdiction. The court chosen is the Court of Modena (Italy).

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**SORIN GROUP USA, Inc**