SORIN | SYNERGY

PH.I.S.I.O.


011655/00 - 02/2009

GB – ENGLISH - DIRECTION FOR USE	5	NL – NEDERLANDS - GEBRUIKSAANWIJZINGEN.....	47
IT – ITALIANO – ISTRUZIONI PER L'USO	11	SE - SVENSKA - BRUKSANVISNING.....	53
FR – FRANCAIS - MODE D'EMPLOI	17	DK - DANSK - BRUGSANVISNING	59
DE - DEUTSCH - GEBRAUCHSANWEISUNG.....	23	FI - SUOMI - KÄYTTÖOHJEET	65
ES – ESPAÑOL - INSTRUCCIONES DE UTILIZACIÓN ..	29	NO – NORSK – BRUKSANVISNING	71
PT – PORTUGUÊS – INSTRUÇÕES DE UTILIZAÇÃO ...	35	CZ - ČEŠTINA - NÁVOD K POUŽITÍ	77
GR – ΕΛΛΗΝΙΚΑ - ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ.....	41		

SORIN GROUP
AT THE HEART OF MEDICAL TECHNOLOGY

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TABLE OF CONTENTS

- A. DESCRIPTION
- B. TECHNICAL FEATURES
- C. INTENDED USE
- D. SAFETY INFORMATION
- E. PREPARATION AND SETUP
- F. PRIMING AND RECIRCULATION PROCEDURE
- G. STARTING PERFUSION
- H. DURING PERFUSION
- I. TERMINATING PERFUSION
- J. BLOOD RECOVERY AFTER PERFUSION
- K. DEVICE CHANGE-OUT
- L. MEDICAL DEVICES FOR USE WITH SYNERGY
- M. RETURN OF USED PRODUCT
- N. LIMITED WARRANTY

A. DESCRIPTION

SYNERGY is an integrated system for hemodynamic support, gas exchange in extracorporeal perfusion and cardiopulmonary bypass for periods of up to 6 hours. SYNERGY comprises the following elements integrated into one assembly: a microporous hollow-fibre membrane oxygenator and heat exchanger, a venous bubble trap, an arterial filter and a centrifugal pump.

The arterial filter includes a 40 µm screen.

The centrifugal pump is equipped with a vane impeller designed to move blood by centrifugal force.

The venous bubble trap is designed for trapping venous air using a 120 µm screen.

SYNERGY is coated with Phosphorylcholine (Ph.I.S.I.O). Devices coated with Ph.I.S.I.O are used when a coated blood path is desired. The Ph.I.S.I.O coating improves the blood compatibility of the device by reducing platelet adhesion to the coated surfaces.

The device is single-use, non-toxic, non-pyrogenic, supplied STERILE and individually packaged, sterilised by ethylene oxide. The level of ethylene oxide residuals in the device is within the limits established by national regulations in the country of use. The device is available in the following versions: [A] SYNERGY for adult bypass.

B. TECHNICAL FEATURES

Maximum flow rate	8000 ml/min
Membrane type	Microporous polypropylene
Membrane surface area	2.0 m ²
Heat exchanger surface area	0.14 M ²
Arterial filter	
- Filter screen	40 µm
- Type	polyester
- Surface area	0.04 m ²
Bubble Trap	
- Filter screen	120 µm
- Type	polyester
- Surface area	0.01 m ²
Static priming volume (post use recovered)	620 ml
Priming volume (oxygenator module + heat exchanger + arterial filter + centrifugal pump + venous bubble trap)	681 ml
Connections:	
- Venous return	3/8" (9.5 mm)
- Arterial filter outlet	3/8" (9.5 mm)
- Blood cardioplegia outlet	pos-lock
Maximum pressure rating:	
- Blood pathway	750 mmHg (100 kpa/ 1bar / 14,5psi)
- Gas pathway at outlet	0 mmHg (0 kpa/ 0 bar/ 0 psi)
- Water pathway	2250 mmHg (300 kpa/ 3 bar / 43,5 psi)

C. INTENDED USE

SYNERGY is intended to be used in adult surgical procedures requiring extracorporeal gas exchange support and blood temperature control.

SYNERGY must not be used longer than 6 hours. Contact with blood for longer periods is inadvisable.

SYNERGY is intended for use only with the Sorin Group Deutschland

D. SAFETY INFORMATION

Information intended to attract the attention of the user to potentially dangerous situations and to ensure correct and safe use of the device is indicated in the text in the following way:

⚠ WARNING

WARNING indicates serious adverse reactions and potential safety hazards for the practitioner and/or the patient that may occur in proper use or misuse of the device and also the limitations of use and the measures to be adopted in such cases.

⚠ CAUTION

CAUTION indicates any special care to be exercised by a practitioner for safe and effective use of the device.

EXPLANATION OF SYMBOLS USED ON THE LABELS

For single use only (Do not reuse)

Batch code (number)
(reference for product traceability)

Use by (Expiry date)

Date of manufacture

Sterile - Ethylene oxide sterilised

Non Pyrogenic

Contains PHTHALATE

Latex free

Warning: Do not resterilize.

Contents sterile only if package is not opened, damaged or broken

REF

Catalogue (code) number

Attention, see instructions for use

This way up

Fragile; handle with care

Keep away from heat

Keep dry

ea

Quantity

The following is general safety information with the aim of advising the operator in preparing to use the device.

⚠ WARNING

- The User should carefully check the device for leaks during set-up and priming. Do not use if any leak is detected.
- The device must be used in accordance with the instructions for use provided in this manual.
- Also specific information is given in the instructions for use at locations in the text where that information is relevant for correct operation.
- For use by professionally trained personnel only.
- Do not use the device if it has been cracked, dropped or otherwise physically damaged.
- SORIN GROUP ITALIA is not responsible for problems arising from inexperience or improper use.
- FRAGILE, handle with care.
- Keep dry. Store at room temperature.
- Always apply and maintain a correct dose and accurate monitoring of the anticoagulant before, during and after the bypass.
- For single use and for single patient use only: during use the device is in contact with human blood, body fluids, liquids or gases for the purpose of eventual infusion, administration or

introduction into the body and due to its specific design it cannot be fully cleaned and disinfected at the end of use. Therefore, reuse on other patients might cause cross-contamination, infection and sepsis. In addition, reuse increases the probability of product failure (integrity, functionality and clinical effectiveness).

- The device must not undergo any further processing.
- Do not resterilize.
- After use, dispose of the device in accordance with applicable regulations in force in the country of use.
- The device must only be used if STERILE.
- The device contains phthalates. Considering the nature of body contact, the limited contact duration and the number of treatments per each patient; the amount of phthalates which might be released from the device do not raise specific concerns about residual risks. Further information is available on request from Sorin Group Italia.
- Do not use the SYNERGY for a period longer than 6 hours. The system has not been qualified through in vivo, ex vivo, or clinical studies for long term use (longer than 6 hours) as a bridge to transplant, or for pending recovery of the natural heart or for extracorporeal circulation (ECMO). The use of the device may result in pump failure, reduced priming capability, excessive blood trauma, degradation and/or corrosion of blood contacting materials with the consequent possibility of particles passing through the bypass circuit to the patient, leaks, and increased potential for gaseous emboli entering the arterial line.
- SYNERGY is a closed loop, active venous drainage system and does not include a venous reservoir. Care must be taken to strictly monitor the pump rpm, delivery flow rate and circuit/patient pressures such that, in case of massive blood loss or slow blood loss, the pump can be adjusted accordingly.
- This device is intended for one-time, single patient use only. Do not resterilize.
- Do not strike the device, as excessive shock may cause damage to the device and cause device malfunction.
- Ensure that SYNERGY is completely primed and de-aired before proceeding with bypass.
- The SYNERGY is a closed loop, kinetically-assisted drainage system. It is therefore mandatory that the utmost care be taken to avoid the accidental introduction of air into the system.
- High vacuum pressures can be generated in the venous line due to cannulae manipulation/movement, flow obstructions including kinking of the venous line, etc. It is highly recommended that the negative pressure at the inlet of the SYNERGY be monitored in order to avoid situations of excessive negative pressure at the venous inlet of the device. It is imperative that the maximum allowable negative pressure limit of -100 mmHg in the venous line is never exceeded. If the venous inlet pressure exceeds this limit, reduce the pump flow rate, then identify and resolve the problem creating the high negative pressure.
- High vacuum pressures can be generated based on cannulae diameter selection. A minimum venous cannula diameter of 28 French is recommended in order to reduce the risk of generating high vacuum pressures in the venous line.
- Do not operate the pump with the centrifugal pump inlet clamped. A negative pressure would be generated in the pump and air bubbles could form.
- The use of the Sorin Group Deutschland Stöckert Electric Remote Clamp (ERC) and a bubble sensor with SYNERGY is highly recommended.
- The water temperature at the heat exchanger inlet must not exceed 42°C (108 °F).
- The water pressure in the heat exchanger must not exceed 2250 mmHg (300 kPa / 3 bar / 44 psi).
- Do not allow the blood side pressure to exceed 750 mmHg (100 kPa).
- Do not use additives or disinfectants, such as bleach, in the heater/cooler while the SYNERGY is connected to the heater cooler.
- Do not allow the pressure on the blood side of the membrane to fall below the pressure on the gas side of the membrane. Blood compartment pressure must always be maintained higher than gas compartment pressure to prevent formation of gas emboli in the blood compartment.
- SYNERGY should always be placed at a level not higher than the patient.
- Ensure that pump flow rate is always greater than cardioplegia and/or haemoconcentration flow rates. In any event, cardioplegia delivery flow rate must not exceed the cardioplegia device maximum rated flow rate.

- Avoid occluding, even accidentally, the gas outlet line: any over-pressurization of the gas compartment may generate micro-emboli in the blood compartment.
- Check the patency of the gas circuit by circulating gas in the circuit. If the gas circuit is occluded, replace the oxygenator.
- Arrange the surgical cannula in such a way that air is not aspirated into the venous line.
- SYNERGY can only be used with the appropriate holder.
- Do not ventilate the SYNERGY while the arterial and venous lines are clamped, as this can cause evaporation of prime solution across the microporous membrane and result in air in the arterial line.
- Prevent that halogenated liquids such as halothane and fluothane come into contact with the polycarbonate structure of the device. This would cause damage which would compromise the integrity and functionality of the device.
- For further information and/or in case of complaints contact SORIN GROUP ITALIA or the authorised local representative.

CAUTION

- Federal law (U.S.A.) restricts this device to sale by or on the order of a physician.
- Do not operate the pump unprimed. Doing so may damage the pump.
- To aid the prevention of gas or particulate embolization, SORIN GROUP ITALIA recommends the use of safety devices, including bubble sensors and pre-bypass filters, for all procedures using the SYNERGY.
- Entry of massive air into the pump will cause the device to deprime and the blood flow to stop.
- The pump must not be run unattended.
- Inner surfaces of the system are Ph.I.S.I.O. coated; currently SORIN GROUP ITALIA is not aware of any contraindication to use of systems having components treated with Ph.I.S.I.O.

E. PREPARATION AND SET-UP

1) HOLDER SET-UP

Position the holder on a pole and fix it by means of the fast clamp connector. Be sure to orient the holder such that the device will be upright, and check that there is sufficient space for the pump driver.

2) MOUNTING THE SYNERGY ON THE HOLDER

WARNING

- Do not use if the sterile packaging is damaged, unsealed, or has been exposed to moisture or other conditions that would compromise the sterility of the device.
- Check the expiry date on the label. Do not use the device after this date.
- The device must be used immediately after opening the sterile packaging.
- The device must be handled aseptically.

Remove the device from the sterile packaging.

WARNING

- Visually inspect and carefully check the device before use. Transport and/or storage conditions other than those prescribed may have caused damage to the device.
- Particularly, turn the device upside down and verify the integrity of the four fins connecting the arterial filter to the oxygenating module.
- Do not use solvents such as alcohol, ether, acetone, etc., as contact may cause damage to the device.
- Do not allow halogenated liquids such as Halothane and Fluothane to come into contact with the polycarbonate housing of the device. This could cause damage which may compromise the integrity and proper functioning of the device.

Fix the SYNERGY oxygenator onto the holder. Check that the water connector locking lever is in the "OPEN" position.

Make sure that the water connectors of the oxygenator are aligned with the Hansen connectors of the holder.

Insert the oxygenator Hansen connectors and push SYNERGY down into the holder.

Turn the locking lever to the "CLOSED" position.

Slide the pump driver into place and secure it to the pump head of SYNERGY. Verify that the centrifugal pump is correctly engaged with the pump driver locking tab. Close the support collar around the pump driver.

Only when the locking lever shows "CLOSED" the oxygenator is correctly inserted in the holder.

⚠ WARNING

- Manually verify that the arterial filter and oxygenating module are stably connected.

3) HEATER/COOLER SET UP

Connect the water tubes to the holder by means of the female Hansen connectors (SORIN GROUP ITALIA code 09028 or equivalent).

⚠ WARNING

- The use of different connectors from those indicated may cause resistance inside the circuit and reduce the efficiency of the heat exchanger.
- Do not obstruct the hole in the lower cover on the heat exchanger as it is the outlet of the safety channel which helps prevent fluids from crossing from one compartment to another.
- The water temperature at the heat exchanger inlet must not exceed 42°C (108 °F).
- The water pressure in the heat exchanger must not exceed 2250 mmHg (300 kPa / 3 bar / 44 psi).

4) CHECK THE HEAT EXCHANGER

⚠ WARNING

- Check the heat exchanger by circulating water inside the heat exchanger for a few minutes. There should be no leaks from the water compartment or from the safety channel hole. In case of leaks discard the unit.

5) DEVICE SET-UP AND CONNECTIONS

- Once the SYNERGY has been positioned in the holder and the centrifugal pump driver is put in place, complete the setup of the remaining extracorporeal circuit components.
- Verify that the 3/8" (9.5 mm) venous line is connected to the venous inlet port located on the SYNERGY bubble trap (ref.1, fig.1).
- Verify that the 3/8" (9.5 mm) arterial line is connected to the arterial filter outlet port (ref.7, fig.1). Connect the centrifugal pump flow probe to circuit where desired.
- Connect the venous purge line provided with the product to the purge Luer connector at the top of the venous bubble trap (ref.2, fig.2). Verify that a means to aspirate air from the bubble trap purge line (e.g., roller pump) is available during prime and device use.
- Verify that the arterial recirculation/purge line is connected to the recirculation/purge Luer of the arterial filter.
- Place the venous bubble sensor assigned to the Electric Remote Clamp on the tubing line between the SYNERGY Bubble Trap Outlet and Centrifugal Pump Inlet.

⚠ CAUTION

- If oxygenated blood for blood cardioplegia needs to be used, remove the red pos lock and connect the 1/4" (6.4 mm) blood line of the cardioplegia circuit to the arterial blood access port of SYNERGY by means of the adaptor (provided with the product).

The blood cardioplegia outlet is fitted with a safety valve that makes connection of the adaptor possible without leakage or dripping during extracorporeal circulation.

⚠ WARNING

- If using the blood cardioplegia outlet during bypass, the blood line connected to the port must be open and not under pressure in order for fluid to flow into the tubing line.

6) SAMPLING SYSTEM

- Mount the sampling stopcock manifold on the dedicated holder. Verify that the stopcock handles are in the "Off" position.
- Connect the arterial sampling line to the Luer next to the recirculation/purge Luer at the top of the arterial filter.

⚠ WARNING

- The arterial sampling site connector Luer does not include a one-way valve. Verify that there is a one-way valve in the sampling line to prevent accidental air introduction into the arterial line.
- c. Connect the venous sampling line to a venous Luer port.

⚠ CAUTION

- Check the security of all Luer connections. The venous side of the SYNERGY including the venous inlet line, bubble trap and tubing line up to the inlet connector of the

centrifugal pump, is under negative pressure. All accessory lines connected to the venous inlet line of the SYNERGY must be connected tightly in order to prevent accidental introduction of air into the negative pressure side of the device.

7) TEMPERATURE AND PRESSURE MONITORING

- Connect the arterial temperature probe to the arterial temperature probe fitting (ref.14, fig.1), and connect the venous temperature probe to a venous temperature probe fitting.
- To measure the oxygenator inlet pressure, connect a pressure measurement line to the venous Luer close to the heat exchanger inlet.
- To measure oxygenator outlet pressure, connect a pressure measurement line to the Luer port close to the arterial blood access port.
- To measure venous inlet pressure in the venous inlet line, attach a pressure measurement line to a venous Luer port.

8) VENTILATING GAS SUPPLY

Connect the 1/4" (6.4 mm) gas supply line to the gas inlet of the oxygenator.

⚠ CAUTION

- The "gas escape" system was designed to prevent any potential risk of occlusion of the gas outlet; such occlusion would cause immediate passage of air to the blood compartment.
- Do not for any reason occlude the external access holes of the gas compartment (ref. 4).
- SORIN GROUP ITALIA recommends the use of a bubble trap or filter on the arterial line to reduce the risk of emboli transmission to the patient

9) VAPOROUS ANAESTHETICS

If vaporous anaesthetics are used, some method of scavenging this gas should be considered. An integral 3/8" (9.5 mm) gas scavenging fitting is available for this purpose. It is located in the centre of the gas escape system. Connect a piece of 3/8 x 3/32" (9.5 x 2.4 mm) tubing to this fitting and to a negative pressure source with a flow rate in excess of the ventilating gas flow rate.

⚠ WARNING

- The only volatile anaesthetics suitable for this use are isoflurane and sevoflurane.
- The methods adopted for vaporous anaesthetic gas scavenging should not in any way increase or reduce the pressure level at the oxygenator fibres.

F. PRIMING AND RECIRCULATION PROCEDURE

⚠ WARNING

- Do not use alcohol-based priming solutions as they can compromise the structural integrity and functionality of the device.

1) CHECK THE HEAT EXCHANGER

Circulate water through the heat exchanger prior to priming the oxygenator. Check the water system for leaks, proper water temperature range, and adequate water flow. If a water to blood leak exists, water will collect in the tubing connected to the oxygenator inlet/centrifugal pump outlet port or drip from the hole in the bottom of the oxygenator module indicating the device must not be used. The maximum pressure rating of the heat exchanger is 2275 mmHg (300 kPa/ 44psi).

2) KEEP THE GAS FLOW CLOSED

3) VERIFY THAT THE RECIRCULATION/PURGE LINE STOPCOCK IS OPEN

4) OCCLUDING THE VENOUS AND ARTERIAL LINES

5) CIRCUIT PRIMING

⚠ WARNING

- The pressure level inside the blood compartment of the oxygenating module may not exceed 750 mmHg (100kPa / 1 bar / 14 psi).
- In order to maintain positive hydrostatic pressure during the priming and recirculation phase, keep the priming bag

above the SYNERGY.

- It is recommended to use a pre-bypass filter to capture any particles that may be present in the circuit or in the priming solution.
- a. While not required, flushing the circuit with carbon dioxide gas will aid in priming.
- b. Prepare to prime with the appropriate amount of priming solution.
- c. Open the venous line attached to the inlet of the SYNERGY bubble trap and slowly allow priming solution into the SYNERGY.
- d. Allow the priming solution to prime the pump, heat exchanger, oxygenator and arterial filter.
- e. Verify that the connection tube between the venous bubble trap and the inlet connector of the pump is completely primed and debubbled.
- f. Prime bubble trap purge line and purge air from bubble trap as needed.
- g. Activate centrifugal pump to recirculate priming solution through recirculation/purge line of the arterial filter.

⚠ WARNING

- During this phase with the recirculation/purge line open and the arterial line clamped, never exceed a pump rate of 1000 rpm.
- h. Remove the clamp from the arterial line, increase the flow rate and recirculate the priming fluid through the A/V loop in order to prime and debubble the circuit.
- i. Set the pump flow rate to 6 l/min. Continue to recirculate for 3-5 minutes. While recirculating, check the system for air bubbles and tap the entire circuit to facilitate the removal of air. Purge any air accumulating at the bubble trap top by aspirating from the purge line.
- j. Ensure that the sampling manifold system is primed. The sample system will prime itself during recirculation when the stopcocks are positioned to allow flow through the sample lines and manifold. The handles point to the off position. Flow through the sample system should be stopped when there is not flow in the bypass circuit.
- k. Ensure that all other parts of the circuit are primed, de-aired and appropriately clamped as needed.
- l. Once the circuit and device are primed, stop the centrifugal pump and clamp the arterial and venous lines. Close the recirculation/purge stopcock.
- m. Prior to bypass, clamp the arterial and venous lines, then pass the clamped lines to the surgeon for connection to the respective cannulae.

⚠ CAUTION

- Always apply and maintain a correct dose and accurate monitoring of the anticoagulant before, during and after the bypass.

⚠ WARNING

- During recirculation, ensure a volume source is open to the circuit. Do not purge using the bubble trap purge line unless a volume source is open to the device. Air will be drawn across the membrane and into the device if the SYNERGY is not open to a volume source.
- If the cardioplegia circuit has been connected to the arterial blood access port, check that the circuit has been primed.
- If tubing lines used for pressure measurement are connected, check that the lines have been primed.
- Occlude the arterial outlet line with a clamp some centimetres away from the outlet port.
- Do not apply negative pressure to the arterial blood access port. Negative pressure in the blood compartment may cause gas microemboli formation.

G. STARTING PERFUSION

1) INITIATING FLOW

Start the centrifugal pump.

⚠ WARNING

- Do not operate the pump for extended periods of time if the arterial outlet to the patient is clamped and the purge/recirculation stopcock is closed at the same time. The temperature inside the pump may increase, increasing the risk of cellular damage.
- Do not remove the arterial outlet clamp at the beginning of

extracorporeal circulation until an adequate outlet pressure has been reached in order to prevent backflow. The pump must be operated to provide a pressure higher than the patient systemic pressure and the circuit head pressure. Monitor the pump flow rate, its rpm and the systemic pressure as indices of potential backflow.

- Carefully monitor for signs of occlusion throughout the circuit, paying particular attention to the venous line and centrifugal pump inlet.
- Do not remove the arterial outlet clamp at the beginning of extracorporeal circulation until an adequate outlet pressure has been reached in order to prevent backflow. The pump must be operated to provide a pressure higher than the patient systemic pressure and the circuit head pressure. Monitor the pump flow rate, its rpm and the systemic pressure as indices of potential backflow.

2) STARTING PERFUSION

⚠ WARNING

- Control the flow rate by adjusting the pump speed. Attempting to adjust the blood flow rate by partial clamping of the device outlet may augment the risk of cellular damage.
- When stopping circulation, clamp the device outlet and stop pump rotation immediately, or reduce the pump speed to less than 1000 rpm if the arterial filter recirculation line is open.

Remove clamp from venous line, remove clamp from arterial line, and initiate bypass according to given standard clinical procedure.

3) ESTABLISH GAS FLOW

⚠ WARNING

- Always open the gas flow after the blood flow. The gas/blood flow ratio must never exceed 2:1.

4) MANAGING THE VENOUS BUBBLE TRAP

As described in Paragraph A (Description), the SYNERGY system is equipped with a venous bubble trap. The bubble trap is fitted with a purge line that is manually controlled by the user using a roller pump or other means suitable for removing the air from the bubble trap. Carefully check for air in the bubble trap, which must be evacuated by aspirating it via the purge line.

⚠ CAUTION

- During circulation, the venous bubble trap is under negative pressure. Any air accumulating at the bubble trap top can be removed only by aspirating it via the purge line. A one-way safety valve is provided in the purge line in order to prevent air from flowing back into the system from the purge line.

5) FUNCTIONING OF HEAT EXCHANGER

Check the temperature of the arterial blood to verify the heat exchanger is working.

6) MANAGING THE ARTERIAL FILTER

Evacuate air, if necessary, by positioning the recirculation/purge stopcock in the purge position. Close the stopcock when purge is finished.

7) BLOOD GAS ANALYSIS

After initiation of bypass, the blood gas parameters should be checked via venous and arterial sampling. After sampling, always move the stopcock back into the "closed" position.

⚠ WARNING

- Venous sampling site is under negative pressure. Always close stopcock before removing syringe.
- Arterial sampling site is under pressure. Always insert a syringe into the sampling site before opening the stopcock.

H. DURING PERFUSION

1) CHECKING THE VENOUS DRAINAGE

The venous drainage is a function of the flow rate of the centrifugal pump and is directly controlled via the centrifugal pump.

⚠ WARNING

- If it is necessary to administer anticoagulant to the patient, refer to § 3 "ADMINISTRATION OF MEDICATION/ADDITIVES".
- The difference between water and blood temperature in the heat exchanger should not exceed 15°C. Higher values may result in the formation of microbubbles from

the dissolved blood gases.

- Using a blood collecting bag in some particular and uncommon conditions (high blood volume, bag almost full, continuous shrinking and swelling of the bag) a collapse of the filtering net with partial occlusion of the outlet can occur. This may cause a momentary block of flow from the outlet connector of the bag. In this case a manual pulling of the bag and a manoeuvre to keep the connector in a straight position is recommended.

2) VENOUS AND ARTERIAL SAMPLING

Use the stopcock manifold for venous and arterial sampling as needed. Open the sample system to permit continuous flushing with arterial blood. This eliminates the necessity of using a flush syringe when drawing an arterial sample. Before taking a venous sample, shut off the arterial side of the sampling system and withdraw at least 10 ml from the centre stopcock.

⚠ CAUTION

- Samples must be taken only when the pump is running and there is blood flow through the SYNERGY. If not, the pressure in the blood compartment may decrease, resulting in the possible formation of air bubbles.

⚠ WARNING

- Venous sampling site is under negative pressure. Always close stopcock before removing syringe.
- Arterial sampling site is under pressure. Always insert a syringe into the sampling site before opening the stopcock.

3) ADMINISTRATION OF MEDICATION/ADDITIVES

While on bypass add all small volume drugs to the venous blood through the sample system to ensure thorough mixing.

⚠ WARNING

- Medication/additives must be injected into the device at the inlet of the venous bubble trap in order to trap and remove any air introduced into the SYNERGY during medication/additive injection.

4) PURGING THE BUBBLE TRAP

The bubble trap must be visually checked periodically for the presence of entrained air. If air is observed or detected by a venous air sensor, manually purge the air from the bubble trap.

⚠ WARNING

- If air is detected, reduce the flow rate, isolate and resolve the cause or source of the air (obstruction of venous inlet line, venous cannula occlusion etc.), and purge the air from the bubble trap as soon as possible.

I. TERMINATING PERFUSION

It should be terminated depending on the conditions of each individual patient. Act as follows:

- 1) Turn the gas flow off.
- 2) Slowly reduce the centrifugal pump speed while gradually occluding the arterial line.
- 3) Open the recirculation/purge line to the recirculation position.
- 4) Occlude the venous line.
- 5) Verify a volume source is open to the device.
- 6) Recirculate as required at 1000 rpm maximum.

⚠ CAUTION

- Do not turn the thermocirculator off during the recirculation phase.
- Verify that the cardioplegia circuit connected to the arterial blood access port is properly occluded.

J. BLOOD RECOVERY AFTER PERFUSION

If desired, maximum blood volume return to the patient may be accomplished by using priming fluid in the priming bag when the blood reaches minimum volume. Pump slowly through the device, delivering blood from the device to the patient while simultaneously verifying that the priming bag is not completely emptied.

K. DEVICE CHANGE-OUT

A spare oxygenator/pump must always be available during perfusion. After 6 hours of use with blood or if situations arise such that, on judgment of the person responsible for perfusion, the safety conditions for the patient are compromised (insufficient oxygenator performance, leaks, anomalous blood parameters etc.), proceed as follows to replace SYNERGY.

- 1) Stop the gas flow, turn off the heater/cooler, and occlude the water

lines.

- 2) Place two clamps on the arterial line (5 cm apart) near the integrated filter, and immediately stop the centrifugal pump.
- 3) Place two clamps on the venous line (5 cm apart) near the inlet of the venous bubble trap.
- 4) Disconnect the water lines, the gas line, the arterial and venous sampling lines, the recirculation/purge line of the integrated arterial filter and the purge line of the venous bubble trap.
- 5) Cut the arterial and venous lines between the clamps.
- 6) Remove the SYNERGY from the holder and disconnect the centrifugal pump from the driver.
- 7) Prepare the SYNERGY unit to be used as replacement, disconnecting the arterial and venous lines using aseptic techniques in order to maintain device sterility.
- 8) Place the new SYNERGY on the holder and follow the above described procedure to mount the device (Ref. § E 2). Connect the pump driver after positioning the replacement SYNERGY properly. Connect the water lines and check integrity of the heat exchanger (Ref. § E 3 and E4).
- 9) Connect the arterial and the venous lines.
- 10) Reconnect the gas line, purge/recirculation lines of the integrated arterial filter, the purge line of the venous bubble trap and the arterial/venous sampling lines.
- 11) Make sure that the purge/recirculation stopcock handle is in the recirculation position.
- 12) Prime the device per § F5 using the AV bypass loop.
- 13) Ensure that all the air is evacuated from the SYNERGY arterial outlet. Close the purge/recirculation line.
- 14) When debubbling has been completed, remove the clamps and resume use of the system.

L. MEDICAL DEVICES FOR USE WITH SYNERGY

The centrifugal pump is intended for use only with Sorin Group Deutschland Stöckert Centrifugal Pump Consoles. Refer to the console operator's manual for console operating procedures.

The Sorin Group Deutschland Stöckert Electric Remote Clamp (ERC) is intended for use only with Sorin Group Deutschland Stöckert Centrifugal Pump Consoles. Refer to the ERC operator's manual for ERC operating procedures.

SYNERGY Holder code 050531.

Stopcock holder

The circuit connections must be made with tubes of a diameter compatible with the dimensions of the connectors located on the device (3/8", 1/4", 1/2").

Temperature control must be carried out with SORIN GROUP ITALIA probes, code 050122 or YSI Series 400 compatible.

Use a Sechrist air/oxygen blender (Sorin Group Italia code 09046) or a system with compatible technical features to control the blood gas concentration.

Any heating/cooling system (thermocirculator) may be used, provided that the connectors to the oxygenator are the Hansen type (Sorin Group Italia code 09028).

Currently SORIN GROUP ITALIA is not aware of any contraindications to use of the device with occlusive or non-occlusive pumps. The use of other types of pumps must be agreed with SORIN GROUP ITALIA.

M. RETURN OF USED PRODUCT

Should the user be dissatisfied with anything related to the quality of the product, the product distributor or the authorized local SORIN GROUP ITALIA representative should be notified.

All parameters considered critical by the user must be reported with particular care and urgency. The following is the minimum information that should be provided:

- Detailed description of the event and, if pertinent, the conditions of the patient;
- Identification of the product involved;
- Lot number of the product involved;
- Availability of the product involved;
- All the indications the user considers useful in order to understand the origin of the elements of dissatisfaction.

Sorin Group Italia reserves the right to authorize, if necessary, recall of the product involved in the notification for assessment. If the product to be returned is contaminated, it must be treated, packed and handled in conformity with the provisions of the legislation in force in the country where the product was used.

⚠ CAUTION

It is the responsibility of the health care institution to adequately prepare and identify the product for return shipment. Do not return products that have been exposed to blood-borne infectious diseases.

N. LIMITED WARRANTY

This Limited Warranty is in addition to any statutory rights of the Purchaser pursuant to applicable law.

SORIN GROUP ITALIA warrants that all reasonable care has been taken in the manufacture of this medical device, as required by the nature of the device and the use for which the device is intended.

SORIN GROUP ITALIA warrants that the medical device is capable of functioning as indicated in the current instructions for use when used in accordance with them by a qualified user and before any expiry date indicated on the packaging.

However, SORIN GROUP ITALIA cannot guarantee that the user will use the device correctly, nor that the incorrect diagnosis or therapy and/or that the particular physical and biological characteristics of an individual patient, do not affect the performance and effectiveness of the device with damaging consequences for the patient, even though the specified instructions for use have been respected.

SORIN GROUP ITALIA, whilst emphasizing the need to adhere strictly to the instructions for use and to adopt all the precautions necessary for the correct use of the device, cannot assume any responsibility for any loss, damage, expense, incidents or consequences arising directly or indirectly from the improper use of this device.

SORIN GROUP ITALIA undertakes to replace the medical device in the event that it is defective at the time of placing on the market or whilst being shipped by SORIN GROUP ITALIA up to the time of delivery to the final user unless such defect has been caused by mishandling by the purchaser.

The above replaces all other warranties explicit or implicit, written or verbal, including warranties of merchantability and fitness for purpose. No person, including any representative, agent, dealer, distributor or intermediary of SORIN GROUP ITALIA or any other industrial or commercial organization is authorized to make any representation or warranty concerning this medical device except as expressly stated herein. SORIN GROUP ITALIA disclaims any warranty of merchantability and any warranty of fitness for purpose with regard to this product other than what is expressly stated herein. The purchaser undertakes to comply with the terms of this Limited Warranty and in particular agrees, in the event of a dispute or litigation with SORIN GROUP ITALIA, not to make claims based on alleged or proven changes or alterations made to this Limited Warranty by any representative, agent, dealer, distributor or other intermediary.

The existing relations between the parties to the contract (also in the case that it is not drawn up in writing) to whom this Warranty is given as well as every dispute related to it or in any way connected to it as well as anything related to it or any dispute concerning this Warranty, its interpretation and execution, nothing excluded and/or reserved, are regulated exclusively by the Italian law and jurisdiction. The court chosen is the Court of Modena (Italy).

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SORIN GROUP USA, Inc