
SORIN

XTRA

GB – ENGLISH	XRES BLOOD COLLECTION RESERVOIR AND COLLECTION SET - INSTRUCTIONS FOR USE.....	3
IT - ITALIANO	XRES BLOOD COLLECTION RESERVOIR E COLLECTION SET - ISTRUZIONI PER L'USO.....	8
FR – FRANÇAIS	XRES BLOOD COLLECTION RESERVOIR ET COLLECTION SET - INSTRUCTIONS D'UTILISATION.....	13
DE - DEUTSCH	XRES BLOOD COLLECTION RESERVOIR UND COLLECTION SET- GEBRAUCHSANWEISUNG.....	18
ES – ESPAÑOL	XRES BLOOD COLLECTION RESERVOIR Y COLLECTION SET - INSTRUCCIONES DE USO.....	23
NL - NEDERLANDS	XRES BLOOD COLLECTION RESERVOIR EN COLLECTION SET - GEBRUIKSINSTRUCTIES.....	28
FI - SUOMI	XRES T BLOOD COLLECTION RESERVOIR JA COLLECTION SET TX - KÄYTTÖOHJE.....	33

CATALOGUE N.	PRODUCT DESIGNATION	CONTENT	PRODUCT DESCRIPTION
04258	XRES T Blood Collection Reservoir	Blood Collection Reservoir TOP	HARD SHELL CONTAINER EQUIPPED WITH AN INTEGRATED FILTER, TO COLLECT BLOOD RECOVERED FROM THE FIELD REMOVING AGGREGATES LARGER THAN 40 µm DIAMETER. IN THE TOP OUTLET VERSION, THE OUTLET PORT IS LOCATED ON THE LID OF THE DEVICE
04259	XRES B BLOOD COLLECTION RESERVOIR	Blood Collection Reservoir BOTTOM	HARD SHELL CONTAINER EQUIPPED WITH AN INTEGRATED FILTER, TO COLLECT BLOOD RECOVERED FROM THE FIELD REMOVING AGGREGATES LARGER THAN 40 µm DIAMETER. IN THE BOTTOM OUTLET VERSION, THE OUTLET PORT IS LOCATED ON THE BOTTOM OF THE DEVICE
04260	Collection Set TX	Blood Collection Reservoir TOP Aspiration Line Vacuum Line	BLOOD COLLECTION RESERVOIR TOP PACKED TOGETHER WITH THE ACCESSORY LINES REQUIRED TO RECOVER BLOOD FROM THE OPERATING FIELD. THE ASPIRATION LINE ALLOWS TO TRANSFER THE BLOOD ASPIRATED FROM THE OPERATING FIELD TO THE RESERVOIR AND CONTEMPORARILY ANTICOAGULATE IT THE VACUUM LINE ALLOWS TO CONNECT THE RESERVOIR TO THE VACUUM SOURCE
04265	Collection Set BX	Blood Collection Reservoir BOTTOM Aspiration Line Vacuum Line	BLOOD COLLECTION RESERVOIR BOTTOM PACKED TOGETHER WITH THE ACCESSORY LINES REQUIRED TO RECOVER BLOOD FROM THE OPERATING FIELD. FOR ASPIRATION AND VACUUM LINES SEE ABOVE FOR DESCRIPTION.
04266	Collection Set TX Cardio	Blood Collection Reservoir TOP Aspiration Line Vacuum Line Cardio Kit	SAME AS THE COLLECTION SET TOP, PLUS AN ADDITIONAL ACCESSORY KIT FOR CONNECTION TO AN OXYGENATOR OR CIRCUIT FOR ECC (CARDIO KIT) THE CARDIO KIT CONTAINS: 1. Y ADAPTER 2. ¼ " ADAPTER TUBING 3. OXYGENATOR EXTENSION LINE 4. LUER-LOCK ADAPTER
04267	Collection Set BX Cardio	Blood Collection Reservoir BOTTOM Aspiration Line Vacuum Line Cardio Kit	SAME AS THE COLLECTION SET BOTTOM, PLUS AN ADDITIONAL ACCESSORY KIT FOR CONNECTION TO AN OXYGENATOR OR CIRCUIT FOR ECC (CARDIO KIT) FOR CARDIO KIT SEE ABOVE FOR DESCRIPTION.

All these devices are single use, non toxic, non pyrogenic, supplied **STERILE** and packaged singly. Sterilized by ethylene oxide. The level of ethylene oxide residuals in the device is within the limits established by national regulations in the country of use.

INTENDED USE

The intended use for the XRES Blood Collection Reservoir, Collection Set Cardio and Collection Set are respectively:

XRES Blood Collection Reservoir

The XRES Blood Collection Reservoir is intended for intra and postoperative sterile recovery and filtration of autologous blood.

Collection Set Cardio

The Collection Set Cardio is intended for intra and postoperative sterile recovery and filtration of the blood shed in **cardiovascular** operations. It is to be used in combination with the aspiration line and with the vacuum line that are included in the set.

The Collection Set Cardio is provided with the Cardio Kit, which includes a set of adapters for connection of the extra-corporeal circulation (ECC) circuit to the Bowl Set of the XTRA Autotransfusion System or to the XRES Blood Collection Reservoir. The

adapters contained in the Cardio Kit allow to extend the use of the autotransfusion techniques to the following procedures:

Concentration of priming fluid during ECC

Recovery of priming fluid residue from the ECC circuit

Collection Set

The Collection Set is intended for intra and postoperative sterile recovery and filtration of the blood shed in **non-cardiovascular** operations. It is to be used in combination with the AAL Aspiration and Anticoagulation Line and with the VEL Vacuum Extension Line that are included in the set.

The products listed above should be used in combination with medical devices listed in section "Medical devices to be used with the products".

CONTRAINDICATIONS

There are no known contraindications for this product. However, the use of blood from this device for reinfusion may be contraindicated, e.g., in the presence of certain potentially harmful substances (refer to cell separator manual for more detail). The responsibility for the use of the device in all cases belongs solely to the physician in charge.

SAFETY INFORMATION

Information intended to attract the attention of the user to potentially dangerous situations and to ensure correct and safe use of the device is indicated in the text in the following way:

⚠ WARNING

Indicates serious adverse reactions and potential safety hazards for practitioner and/or patient that can occur in the proper use or misuse of the device, and also the limitations of use and the measures to be adopted in such cases.

⚠ CAUTION

Indicates any special care to be exercised by a practitioner for the safe and effective use of the device.

EXPLANATION OF THE SYMBOLS USED ON THE LABELS

For single use only (Do not reuse)

Batch code (number)
(reference for product traceability)

Use by (Expiry date)

Date of manufacture

Manufactured by

Sterile - Ethylene oxide sterilised

Non Pyrogenic

Contains PHTHALATE

Latex free

CAUTION: Federal law (usa) restrict this device to sale by or on the order of a physician

Warning: Do not resterilize.

Contents sterile only if package is not opened, damaged or broken

Catalogue (code) number

Attention, see instruction for use

This way up

Fragile; handle with care

Keep away from heat

Keep dry

ea Quantity

The following is general safety information with the aim of advising the operator in preparing to use the device. Also, specific safety information is given in the instructions for use at locations in the text where that information is relevant for correct operation.

⚠ WARNING

The device must be used in accordance with the instructions for use provided in this manual.

The device is intended to be used by professionally trained personnel.

Do not use the device if it has been cracked, dropped or otherwise physically damaged.

SORIN GROUP ITALIA is not responsible for problems arising from inexperience or improper use.

FRAGILE, handle with care.

Avoid any condition which may cause the blood temperature to exceed 37°C (98.6°F).

Keep dry. Store at room temperature.

Always apply and maintain a correct dose and accurate monitoring of the anticoagulant.

For single use and for single patient use only: during use the device is in contact with human blood, body fluids, liquids or gases for the purpose of eventual infusion, administration or introduction into the body and due to

its specific design it cannot be fully cleaned and disinfected at the end of use. Therefore the reuse on other patients might cause cross-contamination, infection and sepsis. In addition, the reuse increases the probability of product failure (integrity, functionality and clinical effectiveness).

The device must not undergo any further processing.

Do not resterilize.

After use, dispose of the device in accordance with applicable regulations in force in the country of use.

The device must only be used if STERILE.

The device contains phthalates. Considering the nature of body contact, the limited contact duration and the number of treatments per each patient; the amount of phthalates which might be released from the device do not raise specific concerns about residual risks. Further information is available on request from Sorin Group Italia.

For further information and/or in case of complaint contact SORIN GROUP ITALIA or the authorized local representative.

⚠ CAUTION

For a detailed description of the circuits, refer to the user's manual of the SORIN GROUP ITALIA cell separator

Refer to the cell separator user's manual for complete instructions on use of the device

Federal law (U.S.A.) restricts this device to sale by or on the order of a physician

SET UP

A XRES Blood Collection Reservoir can be mounted on the integrated holder present on XTRA Autotransfusion System or on an IV pole

HOLDER SET UP

Integrated reservoir holder

Swing open the reservoir holder and raise it to the desired height.

DISPOSABLE SET UP

⚠ WARNING

Sterility is guaranteed only if the sterile packaging is not wet, opened, damaged or broken. Do not use the device if sterility cannot be guaranteed.

Check the expiry date on the label attached. Do not use the device after the date shown.

The device must be used immediately after opening the sterile packaging.

The device must be handled aseptically.

Carry out a visual inspection and carefully check the device before use. Transport and/or storage conditions other than those prescribed may have caused damage to the device.

Do not use solvents such as alcohol, ether, acetone, etc. to prevent damage to the device.

Make the connections using a suitable aseptic technique.

Due to the possibility of operator exposure to bloodborne pathogens (such as HIV, hepatitis viruses, bacteria, Cytomegalovirus, etc.) when handling extracorporeal blood circuits, adequate precautions should be taken at all times to prevent the exposure to and transmission of such agents.

Carefully observe this set for leaks before and during use. Leakage may result in loss of sterility or loss of blood and/or fluid. If leakage is observed before or during use, replace leaking component or retighten leaking connection, as appropriate.

Set up of the XRES Blood Collection Reservoir

The setup of the Blood Collection Reservoir can be performed inserting it into the integrated holder.

1. Remove the reservoir from its package.
2. Keep the locking clip on the reservoir holder open, slide the edge of the reservoir lid into the slot of the holder and push the reservoir into position. Release the locking clip.
3. Tighten all the caps and luer locks on the reservoir lid. In case you are using a Blood Collection Reservoir BOTTOM, close the clamp on the bottom outlet port.

⚠ WARNING

Before use, check that the blood collection reservoir is completely inserted in its housing on the holder. In the event of incorrect installation, any accidental knocks against the reservoir might cause it to detach from the holder

Set up of the AAL Aspiration and Anticoagulation Line

1. Raise the IV pole completely.
2. Prepare a bag or bottle of anticoagulant solution containing heparin, ACD-A or CPDA-1. If heparin is used, mix 30.000-60.000 I.U. of heparin per litre of sterile (injectable) isotonic saline to guarantee sufficient anticoagulation of recovered blood.
3. Remove the aspiration and anticoagulation line from its package and carefully open the outer wrap.
4. Pass the inner sterile wrapped assembly to the sterile field.

- At the sterile field, unwrap the assembly, remove the protective cover, and attach the connector to a suction cannula.
- Pass the other end (the split end) of the aspiration line back from the sterile field to the XTRA Autotransfusion System operator.
- Close the roller clamp on the small bore tubing.
- Remove the protective cover on the large bore tubing of the aspiration and anticoagulation line.
- Attach the tubing to one of the inlet ports on the lid of the reservoir. The inlet ports are 45° angled, orange-capped, marked SUCTION.
- Hang the anticoagulant solution container on one of the holders of the IV pole.
- Using aseptic technique, spike the container of the anticoagulant solution.

Set up of the VEL Vacuum Extension Line

- Attach one end of the vacuum line to the XVAC Vacuum System (or another vacuum source)
- Attach the other end to the vacuum port (yellow capped, marked VACUUM) on the lid of the XRES Blood Collection Reservoir.

Set up of the Cardio Kit (for Collection Set Cardio only, otherwise skip next steps)

Connecting the Y adapter to the XRES Blood Collection Reservoir

- Close the clamps on the Y adapter branches waiting for subsequent connections.

NOTE: The branch ending with a male port is used to connect the reservoir to the Bowl Set, the branch ending with a female port is used to connect the XRES Blood Collection reservoir to a circuit or oxygenator for ECC through the Oxygenator Extension Line.

- Connect the Y adapter ¼" port to the reservoir outlet port:
 - located on the reservoir lid, marked WASH and covered with a blue cap on the XRES T Blood Collection Reservoir
 - located on the bottom of the XRES B Blood Collection Reservoir

Connecting the Oxygenator Extension Line to a circuit or oxygenator for ECC

- Before adding the priming liquid to the ECC circuit, connect the 1/4" Adapter Tubing contained in the Cardio Kit to a 1/4" connector of the circuit or oxygenator for ECC.
- Close the ¼" Adapter Tubing with a clamp, keeping the protective cap in place to maintain sterility. During priming circulation remove the air from the circuit.
- Connect one end of the Oxygenator Extension Line to the free end of the ¼" Adapter Tubing
- Connect the other end to the female connector on the Y Adapter branch or directly to the Bowl Set.

INSTRUCTIONS FOR USE

Blood collection and anticoagulation

- Switch on the vacuum system, following below instructions according to the system used
 - If using XVAC vacuum system linked to XTRA Autotransfusion System:
 - Make sure XVAC is connected to the XTRA and the power switch on the back of the vacuum unit is on ON.
 - Turn on the XTRA Autotransfusion System and wait for system's boot up
 - Regulate the vacuum level from the vacuum display present on the cell separator's screen

NOTE:

- Once system's boot up is finished, the level of vacuum is controlled from the vacuum display present on the cell separator's screen
- If during an emergency XVAC is required while the XTRA OS boots up, the aspiration can be controlled directly from the control panel present on the XVAC vacuum system (refer to chapter 13 of XTRA Operator's Manual for detailed information), until the boot up is finished. Then, it will be controlled only from the XTRA screen.

- If using XVAC vacuum system stand alone:
 - Make sure the power switch on the back of the unit is on ON.
 - Turn on XVAC using the ON/OFF switch present on the control panel
 - Regulate the level of vacuum from the control panel
- If using a vacuum regulator
 - Switch on the vacuum regulator
 - Close all the lines apart from the vacuum line with a clamp (the aspiration line and the Y adapter if present) and regulate the vacuum level.
 - Then, reopen the lines, keeping the vacuum on.

WARNING

SORIN GROUP ITALIA recommends not to exceed a negative pressure of -150 mmHg (20 Kpa/0,2 Bar/-2,9 psi) during recovery procedure as an excess vacuum increases hemolysis of the recovered blood

- Regulate the vacuum at a level no greater than 150 mmHg (20 kPa) in accordance with the AABB guidelines.
- Verify that the protective cover on the end of the suction line in the sterile field has been removed and that the aspiration tip is open to the atmosphere (not blocked).
- With the vacuum system ON, open the roller clamp on the small bore tubing and allow 200 ml of anticoagulant solution to be drawn into the reservoir before starting blood collection. This volume should ensure adequate wetting of the blood contact surfaces. Close the roller clamp if there will be any delay before processing.
- After the transfer of the anticoagulant solution to the reservoir, adjust the flow rate of the anticoagulant solution to between 60 and 100 drops a minute. The lower flow rate is suitable for ACD-A or CPDA-1 and the higher rate for heparin. In case of heavy bleeding, appropriately increase the rate.

NOTE: The flow of anticoagulant into the reservoir is manually controlled by the roller clamp on the anticoagulant line. This flow must be adjusted to the rate of blood collection from the surgical field. If the rate of collection varies without adjusting the flow of anticoagulant, the ratio of anticoagulant to blood can be too low or too high. The blood in the reservoir may clot if there is too little anticoagulant.

WARNING

In the event of excess heparin in collection reservoir due to inappropriate ratios, the salvaged blood may contain residual heparin.

In the event of decreased patient antithrombin III levels if using heparin anti-coagulation, consult the physician in order to provide alternate anticoagulation.

WARNING

When in the need of transferring into the reservoir unfiltered blood coming from another source, use the filtered luer lock connection.

CAUTION

Failure to maintain adequate anticoagulation during blood collection can cause excessive clotting and possible clogging of the collection reservoir or centrifuge bowl.

Prior to collecting any blood, prime the reservoir with approximately 200 mL of the anticoagulant solution.

When an anticoagulant solution using ACD-A or CPDA-1 is used during blood recovery, do not use ringer lactate solutions in subsequent wash.

The XRES Blood Collection Reservoir is equipped with three graduated scales 180° angled to allow the operator to monitor the level of fluid present inside the reservoir from any position. The indication provided by the graduated scales must be considered only as a gross indication of the fluid contained in the device.

The XRES Blood Collection Reservoir is equipped with a Safety level valve device for interrupting aspiration and preventing the blood collected from entering the vacuum line. Should this safety level valve activate, follow the instructions below to reset aspiration to the operative field:

- Temporarily switch off the vacuum and process the fluid contained in the reservoir until the level is below the safety valve.
- Pull the yellow tab located on the connector marked VACUUM and check that the valve has dropped down in the guide.
- Release the yellow tab and reset aspiration.

WARNING

A pressure relief valve is provided to prevent implosion of the device if exposed to a sudden, deep vacuum

- The negative pressure inside the XRES Blood Collection Reservoir must not exceed -300 mmHg (40 KPa/0,4 bar/5,7 psi).
- Check the cap seals of the connectors not in use, pushing them fully on.
- Do not occlude the safety valve release present on the XRES Blood Collection Reservoir with foreign bodies to prevent the potential possibility of implosion

Using the XRES Blood Collection Reservoir connected to a Bowl Set

- When ready to start processing, connect the Bowl Set inlet port (on the fill line, blue color) to the XRES Blood Collection Reservoir outlet port:
 - marked WASH and covered with a blue cap, for a XRES T Blood Collection Reservoir
 - located on the bottom of the XRES T Blood Collection Reservoir
- In case a Y adapter is used (Collection Set Cardio), connect the Bowl Set inlet port (on the fill line, blue color) to the Y adapter male branch and open the corresponding clamp. Keep the clamp on the other branch of the Y adapter (for connection to the Oxygenator Extension Line) closed.
- Proceed with further steps for blood processing according to the instructions for use contained in the operator's manual of the cell separator

Concentration of priming fluid during ECC (for Collection Set Cardio only)

- Isolate the Y adapter from the reservoir, by closing the 1/4" tubing segment connecting the Y adapter to the reservoir with a clamp. Open the clamps located on the Y adapter branches.
- Remove the clamp used on the 1/4" Adapter Tubing previously connected to the ECC circuit and open any other clamp that may be present on the

- Oxygenator Extension Line.
3. Whenever the Oxygenator Extension Line is directly connected to the Bowl Set, ensure that all the clamps are open.
4. Transfer the priming fluid to the cell separator device following the instructions for use of the device, taking care not to create any negative pressure within the ECC circuit.
5. During the return of the concentrated priming fluid, avoid any air being routed to the ECC circuit.

Recovery of priming residue from the ECC circuit (for Collection Set Cardio only)

During cardiovascular operations, it is advisable to transfer the priming residue from the ECC circuit to the reservoir, to minimize the time the pump console remains in the operating theatre.

1. Clamp the aspiration line
2. Open both clamps on the Y adapter branches. Make sure the ¼" tubing segment connecting the Y Adapter to the reservoir is open. Remove the clamp used on the 1/4" Adapter Tubing previously connected to the ECC circuit and open any other clamp that may be present on the Oxygenator Extension Line
3. Transfer the priming residue to the reservoir from the ECC circuit until it is completely empty, appropriately adjusting the vacuum.
4. After transferring the priming residue to the reservoir, close all the clamps between the Y Adapter and the ECC circuit and remove the connection to the circuit. The pump console is now free of connections and the blood processing can start.
5. Reopen the aspiration line

⚠ CAUTION

Wash the recovered blood with at least the volume of saline solution recommended in the instructions for use of SORIN GROUP ITALIA cell separators for autotransfusion, to guarantee removal of excess anticoagulant and other undesired components collected.

REINFUSING SALVAGED BLOOD

⚠ WARNING

It is recommended that intraoperative and postoperative autologous blood be processed. The washing operation is necessary whenever the recovery is carried out in the presence of heavy haemolysis caused by surgery, or in cases of high coagulopathy or hepatic or renal failure.

The decision to re-infuse the blood components obtained by this concentration-washing procedure is the responsibility of the attending physician.

Do not reinfuse the patient's blood from the primary reinfusion bag when it is connected to the cell separator's wash circuit. Reinfusion from the primary reinfusion bag connected to the circuit could lead to air embolism to the patient.

Do not reinfuse under pressure (i.e., do not use a blood pressure cuff on the secondary reinfusion bag). Reinfusion under pressure could lead to air embolism.

To minimize the complications of particulate matter infusion, use of an in-line microaggregate filter on the patient reinfusion line is highly recommended. Failure to use an in-line filter may result in reinfusing particles potentially harmful to the patient.

To reduce risk of air embolism, remove all air from the primary reinfusion bag before handing the bag over for reinfusion.

The American Association of Blood Banks recommends the following guidelines for expiration of salvaged blood¹.

- a. If not transfused immediately, units collected and processed under sterile conditions with a device for intraoperative blood collection that washes with 0.9% saline, USP, shall be stored under one of the following conditions prior to completion of transfusion:

At room temperature for up to 4 hours after completion of processing

at 1 °C to 6 °C for up to 24 hours, provided that storage at 1°C to 6 °C is begun within 4 hours after completion of processing

- b. The transfusion of shed blood collected under postoperative or posttraumatic conditions shall be completed within 6 hours of initiating the collection.

1. American Association of Blood Banks. *Standards for Perioperative Autologous Blood Collection and Administration*. 3rd ed. Bethesda, MD: American Association of Blood Banks, 2007.

Using XRES Blood Collection Reservoir Quick Reinfusion Port

XRES T Blood Collection Reservoir

The XRES T Blood Collection Reservoir is provided with a pierceable port on the bottom, to gravity drain blood out of the reservoir. This port can be accessed with a spike connector and may be used to empty the reservoir in case needed, or to quickly reinfuse the unprocessed salvaged blood to the patient, in emergency situations.

Observe the following instructions to access the Quick Reinfusion Port:

1. Shut off the vacuum source, pull the yellow tab on the vacuum port of the reservoir lid, then raise the reservoir as high as needed
2. Insert a spike through the integral membrane
3. In case of quick reinfusion, completely prime the line before connection to the patient

XRES B Blood Collection Reservoir

The XRES B Blood Collection Reservoir is provided with a ¼" male outlet port on the bottom, to gravity drain blood out of the reservoir. This port may be used to empty the reservoir in case needed, or to quickly reinfuse the unprocessed salvaged blood to the patient, in emergency situations.

Observe the following instructions to use the bottom outlet port as a Quick Reinfusion Port:

1. Shut off the vacuum source, pull the yellow tab on the vacuum port of the reservoir lid, then raise the reservoir as high as needed
2. Connect an infusion line to the ¼" bottom outlet port male connector, using an adapter if needed
3. In case of quick reinfusion, completely prime the line before connection to the patient

⚠ WARNING

Limit the reinfusion of unprocessed blood to the surgical situations of real necessity. The decision and responsibility of reinfusing unprocessed blood belongs solely to the physician in charge. When reinfusing unprocessed blood, the use of a microaggregate filter on the patient reinfusion line is mandatory.

MEDICAL DEVICES TO BE USED WITH THE PRODUCTS

The XRES Blood Collection Reservoir, Collection Set Cardio and Collection Set must be connected to the Sorin group italia Bowl Sets for removal of undesired substances collected from the surgical field.

The connections to the product must be of a diameter which is compatible with the dimensions of the connectors on the device (3/8" or 1/4").

The aspiration must be performed with SORIN GROUP ITALIA surgical aspirator or other aspirators having compatible technical characteristics.

⚠ CAUTION

The User should observe the warnings and cautions and follow Instruction for Use accompanying the separate device.

RETURN OF USED PRODUCTS

Should the user be dissatisfied with anything related to the quality of the product, the product distributor or the authorized local SORIN GROUP ITALIA representative should be notified.

All parameters considered critical by the user must be reported with particular care and urgency. The following is the minimum information that should be provided:

- Detailed description of the event and, if pertinent, the conditions of the patient;
- Identification of the product involved;
- Lot number of the product involved;
- Availability of the product involved;
- All the indications the user considers useful in order to understand the origin of the elements of dissatisfaction.

SORIN GROUP ITALIA reserves the right to authorize, if necessary, recall of the product involved in the notification for assessment. If the product to be returned is contaminated, it must be treated, packed and handled in conformity with the provisions of the legislation in force in the country where the product was used.

⚠ CAUTION

It is the responsibility of the health care institution to adequately prepare and identify the product for return shipment. Do not return products that have been exposed to blood borne infectious diseases.

ONLY for US customers

If for any reason the product must be returned to the manufacturer, a returned good authorisation (RGA) number is required from Sorin Group USA, Inc. prior to shipping. If the product has been in contact with blood or blood fluids, it must be thoroughly cleaned and disinfected before packing. It should be shipped in either the original carton or an equivalent carton to prevent damage during shipment, and it should be properly labelled with an RGA number and an indication of the biohazardous nature of the content in the shipment.

The shipping address for returned goods in the US is:

SORIN GROUP USA, Inc.
Returned CV Products
14401 West 65th Way
Arvada, CO 80004-3599
FAX (800) 323 4031

LIMITED WARRANTY

This Limited Warranty is in addition to any statutory rights of the

Purchaser pursuant to applicable law.

SORIN GROUP ITALIA warrants that all reasonable care has been taken in the manufacture of this medical device, as required by the nature of the device and the use for which the device is intended. SORIN GROUP ITALIA warrants that the medical device is capable of functioning as indicated in the current instructions for use when used in accordance with them by a qualified user and before any expiry date indicated on the packaging. However, SORIN GROUP ITALIA cannot guarantee that the user will use the device correctly, nor that the incorrect diagnosis or therapy and/or that the particular physical and biological characteristics of an individual patient, do not affect the performance and effectiveness of the device with damaging consequences for the patient, even though the specified instructions for use have been respected. SORIN GROUP ITALIA, whilst emphasizing the need to adhere strictly to the instructions for use and to adopt all the precautions necessary for the correct use of the device, cannot assume any responsibility for any loss, damage, expense, incidents or consequences arising directly or indirectly from the improper use of this device. SORIN GROUP ITALIA undertakes to replace the medical device in the event that it is defective at the time of placing on the market or whilst being shipped by SORIN GROUP ITALIA up to the time of delivery to the final user unless such defect has been caused by mishandling by the purchaser. The above replaces all other warranties explicit or implicit, written or verbal, including warranties of merchantability and fitness for purpose. No person, including any representative, agent, dealer, distributor or intermediary of SORIN GROUP ITALIA or any other industrial or commercial organization is authorized to make any representation or warranty concerning this medical device except as expressly stated herein. SORIN GROUP ITALIA disclaims any warranty of merchantability and any warranty of fitness for purpose with regard to this product other than what is expressly stated herein. The purchaser undertakes to comply with the terms of this Limited Warranty and in particular agrees, in the event of a dispute or litigation with SORIN GROUP ITALIA, not to make claims based on alleged or proven changes or alterations made to this Limited Warranty by any representative, agent, dealer, distributor or other intermediary. The existing relations between the parties to the contract (also in the case that it is not drawn up in writing) to whom this Warranty is given as well as every dispute related to it or in any way connected to it as well as anything related to it or any dispute concerning this Warranty, its interpretation and execution, nothing excluded and/or reserved, are regulated exclusively by the Italian law and jurisdiction. The court chosen is the Court of Modena (Italy).

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