
DIDECO CARDIOPLEGIA DELIVERY SYSTEM

Blood and Clear Cardioplegia Solution Delivery Set

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A. DESCRIPTION

D 512, D 513, D 514, D 516, D 512/S, D 513/S, D 514/S, D 516/S are sets designed to administer cardioplegic solution containing blood with a fixed blood/asanguineous cardioplegic solution ratio.

Cardioplegic solution is used for perfusion of the coronaries during cardiopulmonary by-pass, if the aorta is clamped and the intention is to protect the myocardium during the anoxic phase.

Circuits equipped with selector (/S) allow antegrade and retrograde perfusion of the coronaries.

The sets are disposable, non-toxic and non-pyrogenic and are supplied **STERILE** in individual packs. They are sterilized with ethylene oxide. The residual ethylene oxide level in the sets complies with the requirements of the legislation in force in the country of use.

The sets consist of:

1. Two lines for supply of the asanguineous cardioplegic solution and oxygenated blood which also act as pump segments.
2. D 720/C heat exchanger incorporating bubble chamber and temperature ports.
3. Three-way stopcock with pressure and vent line.
4. Delivery line with male luer-lock connector for connection to coronary access cannula.
5. Two-way selector for antegrade and retrograde perfusion (/S models).

B. TECHNICAL SPECIFICATIONS

The systems differ in the different proportion of blood in relation to cardioplegic solution and whether or not they are equipped with selector.

More specifically:

- D 512 4:1 (4 parts blood and 1 part asanguineous cardioplegic solution); for flow calculation refer to a diameter of 1/4" and 1/8" for blood and asanguineous cardioplegic solution respectively (6.36 and 3.17 mm). Ht obtained 80% of original Ht.
- D 512/S with antegrade/retrograde delivery selector.
- D 513 2:1 (2 parts blood and 1 part asanguineous cardioplegic solution); for flow calculation refer to a

- diameter of 1/4" and 3/16" for blood and asanguineous cardioplegic solution respectively (6.36 and 4.8 mm). Ht obtained 65.5% of original Ht.
- D 513/S with antegrade/retrograde delivery selector.
- D 514 1:1 (1 part blood and 1 part asanguineous cardioplegic solution); for flow calculation refer to a diameter of 1/4" for both tubes (6.36 mm). Ht obtained 50% of original Ht.
- D 514/S with antegrade/retrograde delivery selector.
- D 516 9:1 (9 parts blood and 1 part asanguineous cardioplegic solution); for flow calculation refer to a diameter of 1/4" and 1/12" for blood and asanguineous cardioplegic solution respectively (6.36 and 2.11 mm). Ht obtained 90% of original Ht.
- D 516/S with antegrade/retrograde delivery selector.

C. INTENDED USE

D 512, D 513, D 514, D 516, D 512/S, D 513/S, D 514/S, D 516/S are sets for the administration of cardioplegic solution containing blood, with the aid of a low flow pump. This set allows control of the flow, the temperature and infusion (line) pressure of the cardioplegic fluid.

The pump permits total control of the infusion parameters, flow and quantity of solution delivered.

The perfusion temperature depends on the thermal parameters of the cold water source which supplies the D 720/C.

The formulae and graphs are an aid to extrapolation of the cardioplegic fluid temperature. Sets must not be used for more than 6 hours.

Contact with the perfusion solution for longer periods is not recommended.

Sets must be used in combination with the medical devices listed in point I (Medical devices for use with cardioplegia sets).

D. SAFETY INFORMATION

In the text, information intended to draw the user's attention to the need to prevent hazardous situations and assure correct, safe use of the device are identified as follows:

WARNING

Indicates serious consequences and potential dangers for the user and/or the patient deriving from the use of the device in normal conditions of use or improper use, as well as limitations on its use and the measures to be taken if these events should occur.

CAUTION

Indicates every possible precaution the user must take for safe, effective use of the device.

KEY TO THE SYMBOLS USED ON THE LABELS

For single use only
(Do not reuse)

Batch code (number) (reference for product traceability)

Use by (Expiry date)

Date of manufacture

Sterile - Ethylene oxide sterilised

Non Pyrogenic

Contains PHTHALATE

Latex free

Warning: Do not resterilize.
Contents sterile only if package is not opened, damaged or broken

Catalogue (code) number

Attention, see instruction for use

This way up

Fragile; handle with care

Quantity

Keep away from heat

Keep dry

The following are items of general information about safety, intended to provide the operator of the set with appropriate warnings. In addition, specific information about safety is included in these instructions for use where they may affect the operation to be performed.

WARNING

- The device must be used in accordance with these instructions for use.

- The device is intended for professional use.
- SORIN GROUP ITALIA accepts no responsibility for problems caused by negligence or improper use.
- FRAGILE: handle with care
- Keep dry. Store at room temperature.
- For single use and for single-patient use only. During use the device is in contact with human blood, body fluids, liquids or gases for the purpose of eventual infusion, administration or introduction into the body, and due to its specific design it cannot be fully cleaned and disinfected after use. Therefore, reuse on other patients might cause cross-contamination, infection and sepsis. In addition, the reuse increases the probability of product failure (integrity, functionality and clinical effectiveness).
- The device contains phthalates. Considering the nature of contact with the body, the limited duration of contact and the number of treatments per patient, the amount of phthalates which might be released from the device does not raise specific concerns about residual risks. Further information is available on request from Sorin Group Italia.
- Do not subject to further treatments.
- Do not resterilize.
- After use, dispose of the device in accordance with the applicable regulations in the country of use.
- The device must only be used if STERILE.
- If the device is supplied NON STERILE (in which case the pack will be marked "non sterile"), reach agreement on the sterilisation method with SORIN GROUP ITALIA or its authorised representative.
- It is recommended that a cold agglutinin test be performed with the patient's blood before using the set. The temperature of cardioplegic solution should not be lower than the cold agglutination temperature if the test shows the cold agglutinin to be haemolytic.
- When used in conjunction with membrane oxygenators, to avoid introducing air into the circuit:
 - do not start cardioplegia pump until arterial pump is on
 - do not stop arterial pump until cardioplegia pump is off
 - arterial pump flow must always exceed cardioplegia pump flow.
- For further information and/or in case of complaint, contact SORIN GROUP ITALIA or your local authorised representative.

E. ASSEMBLY **WARNING**

- The cardioplegia sets must be installed before priming the extracorporeal tubing system and the air venting procedure must be carried out before extracorporeal circulation begins.
- Sterility is guaranteed provided the sterile pack

has not been wet, opened, tampered with or damaged.

- Do not use the device if sterility is not guaranteed.
- Check the shelf-life expiry date on the label provided. Do not use the device after this date.
- The device must be used immediately after the sterile pack is opened.
- The device must be handed in an aseptic manner.
- Make a visual check on the device and inspect it carefully before use. Inappropriate transport and/or storage conditions may have damaged the product.
- Do not use solvents such as alcohol, ether, acetone etc.; they may damage the product on contact.
- Do not allow halogenated fluids such as Halothane or Fluothane to come into contact with the device's polycarbonate structure. This would cause damage which might impair the device's operation.

1) FIX THE HEAT EXCHANGER TO THE HOLDER AND CONNECT TO THE THERMOCIRCULATOR

Remove the set from its sterile pack.

Assemble the set following the instructions given on the labels on the ends of the tubes. The D 720/C must be placed in the holder D 611 and connected to the thermocirculator by means of the female Hansen connectors, SORIN GROUP ITALIA code 9028.

2) CHECK THE HEAT EXCHANGER

The D 720/C is checked by circulating water through the exchanger for a few minutes; the structure is guaranteed undamaged if there are absolutely no leaks from the water compartment.

3) OTHER CONNECTIONS

Connect the pressure line to a pressure gauge. Connect the vent line to a sterile bag or to the cardiostomy reservoir of the extracorporeal circuit. Connect two probes mod. 9026 or equivalent model Y.S.I. series 400 compatible probes to the D 720/C temperature ports.

4) INSTALL THE PUMP LINES

Install the blood and cardioplegia lines in the cardioplegia pump head, making sure that the pump segments are fitted in such a way that the flow direction is towards the "Y" fitting.

This pump should be equipped with the correct tube clamping inserts.

WARNING

Ensure that the cardioplegic solution line is on top of the blood line and that both lines are fully occluded.

5) CONNECT THE CARDIOPLEGIC SOLUTION LINE

Connect the cardioplegic solution line to the solution container by spiking the latter.

6) CONNECTING THE BLOOD LINE

Connect the blood line to the coronary outlet (1/4")

of the oxygenator; this connection must always be made before priming the oxygenator.

7) CONNECT THE DELIVERY LINE

Connect the delivery line set, ensuring sterility: the end with the luer-lock connector will be connected to the vascular access prepared (coronary cannulae or needle cannula).

Version with selector

Connect the luer lock of the delivery line to the selector input.

F. PRIMING PROCEDURE

1) PRIMING THE SET

Prime the oxygenator and the extracorporeal tubing system, circulate and debubble before priming the cardioplegia system. Make sure that the pump is completely occlusive and start priming with prime solution from the oxygenator.

Clamp the vent tube, turn the D 720/C upside down and start filling the system by operating the pump at a flow between 20 and 30 ml/min.

2) VENT AIR FROM THE EXCHANGER

As soon as liquid flows from the outlet connector, clamp the outlet tube and return the heat exchanger to its initial position. At the same time, open the vent line on the top of the bubble chamber.

Venting is completed by continuing to pump air from the D 720/C at a flow of 20-30 ml/min while tapping it with a plastic coated object. An air pocket of a few cubic centimetres should be retained to ensure that cardioplegic solution with a high potassium content never reaches the extracorporeal tubing system.

3) VENT AIR FROM THE CIRCUIT

Clamp the D 720/C vent line and at the same time open the delivery line. The section of the circuit which delivers the fluid to the patient will fill and discharge its air; when all air has been vented, clamp this line again.

The surgeon can now make the connection to the vascular access prepared (coronary cannulae or needle cannula) and then vent the air from it. The section of the delivery line immediately upstream of the vascular access should be clamped.

Version with selector

Prime the selector and vent the air.

G. DELIVERY PROCEDURE

WARNING

- Always make sure that the arterial pump is running before operating the cardioplegia pump. This will ensure air is not extracted from the oxygenator.

1) Deliver cardioplegia after removing the clamp from the delivery line and the clamps (if any) from the pump circuit. Switch on the thermocirculator one minute before perfusion begins. Adjust pump speed

to obtain the desired flow and pressure values.

Version with selector

Select antegrade or retrograde perfusion.

Constantly monitor the procedure so that the blood level in the oxygenator venous reservoir and the quantity and quality of the liquid in the asanguineous cardioplegic solution container are always under control.

- 2) To terminate cardioplegia delivery stop cardioplegia pump and clamp delivery line.
- 3) To administer further doses check that there is no air in the circuit and repeat the procedures described in points 1 and 2.

WARNING

- Check carefully that the pressure is within the recommended range; different pressure values may lead to inadequate perfusion or damage to the coronaries.
- The line pressure is not the same as the perfusion pressure since it depends above all on the resistance of the delivery circuit: tube and cannulae. The actual perfusion pressure is generally lower than the reading.
- Keep a constant strict check on the temperature of the cardioplegic solution; this temperature will depend on the overall flow of the asanguineous cardioplegic solution and the temperature of the cooling water.

CAUTION

- Monitor the cardioplegia delivery in terms of quantity, time, pressure and temperature.
- Always make sure that the quantities of fluids for the cardioplegic perfusion procedure are sufficient for the session and to ensure troublefree extracorporeal circulation.
- Make sure that the fluid level in the D 720/C bubble chamber is always close to the maximum value. If this level drops, vent the air by operating the pump at a flow of 20-30 ml/min and conveying the air towards the vent set.
- Replace the asanguineous cardioplegic solution container in good time if the amount it contains appears insufficient for completion of the delivery procedure.
- Always monitor the system closely during the delivery of cardioplegic solution. Air in the circuit constitutes a danger to the patient.

H. REPLACING THE CARDIOPLEGIA SET

A spare set should always be to hand during perfusion. If situations arise in which the person in charge considers the patient is at risk (inadequate performance, leak from hole on D 720 bottom cover, etc.), replace the set as described below.

CAUTION

Use sterile techniques throughout the replacement procedure.

- 1) Stop the pump.
- 2) Place two clamps on the blood line connected to the coronary outlet of the oxygenator (5 cm apart).
- 3) Switch off the thermocirculator, clamp and remove the water lines. Clamp the asanguineous cardioplegic solution line.
- 4) Remove the vent line (D 720 C only) and all the monitoring lines.
- 5) Clamp the D 720/C outlet line distal to the pos-lock.

Version with selector

If the set includes the selector (S sets) clamp the selector inlet line.

- 6) Shut the stopcock on the connector linking the vent and pressure lines and unscrew the pos-lock. These two lines will be used with the new set.
- 7) Cut the blood line at the point between the two clamps and fit a 1/4"-1/4" connector.
- 8) Remove the spike from the asanguineous cardioplegic solution container.
- 9) Remove the D 720/C and the set to be replaced from the holder.
- 10) Fit the new cardioplegia set. Connect all the lines (blood line to the 1/4"-1/4" connector, outlet line, vent line and pressure line). Remove the clamps from the blood line and the asanguineous cardioplegic solution line.

Version with selector

Connect the outlet line to the selector inlet.

- 11) Remove the clamps from the water line and switch on the thermocirculator; check the new D 720/C for leaks.
- 12) Prime the new circuit and vent any air bubbles, as described in the priming procedure.
- 13) Restart perfusion.

I. MEDICAL DEVICES FOR USE WITH THE CARDIOPLEGIA SETS

Temperature monitoring must be by means of a SORIN GROUP ITALIA code 9026 or Y.S.I. series 400 compatible probe.

There are no contraindications for the use of any heating/cooling system (thermocirculator), except that the water connections must be the SORIN GROUP ITALIA code 9028 Hansen type.

To the best of SORIN GROUP ITALIA's current knowledge, there are no contraindications to the use of the device with occlusive or non occlusive roller pumps or centrifugal pumps.

The use of other types of pump must be agreed with SORIN GROUP ITALIA.

To the best of SORIN GROUP ITALIA's current knowledge, there are no contraindications to the use of silicone tubes with occlusive roller pumps for extracorporeal circulation.

Incorrect occlusion setting or incorrect use of the

pump are not covered by the silicone tube quality guarantee, since these are variables which cannot be controlled by SORIN GROUP ITALIA.

J. RETURN OF USED PRODUCTS

Should the user be dissatisfied with anything related to the quality of the product, the product distributor or the authorized local SORIN GROUP ITALIA representative should be notified.

All parameters considered critical by the user must be reported with particular care and urgency. The following is the minimum information that should be provided:

- Detailed description of the event and, if pertinent, the conditions of the patient;
- Identification of the product involved;
- Lot number of the product involved;
- Availability of the product involved;
- All the indications the user considers useful in order to understand the origin of the elements of dissatisfaction.

SORIN GROUP ITALIA reserves the right to authorize, if necessary, recall of the product involved in the notification for assessment. If the product to be returned is contaminated, it must be treated, packed and handled in conformity with the provisions of the legislation in force in the country where the product was used.

CAUTION

It is the responsibility of the health care institution to adequately prepare and identify the product for return shipment. Do not return products that have been exposed to blood borne infectious diseases.

ONLY for US customers

If for any reason the product must be returned to the manufacturer, a returned good authorisation (RGA) number is required from Sorin Group USA cardiovascular Inc. prior to shipping.

If the product has been in contact with blood or blood fluids, it must be thoroughly cleaned and disinfected before packing. It should be shipped in either the original carton or an equivalent carton to prevent damage during shipment, and it should be properly labelled with an RGA number and an indication of the biohazardous nature of the content in the shipment.

The shipping address for returned goods in the US is:
Sorin Group USA, Inc.
Returned CV Products
14401 West 65th Way
Arvada, CO 80004-3599
FAX (800) 323 4031.

K. LIMITED WARRANTY

This Limited Warranty is in addition to any statutory rights of the Purchaser pursuant to applicable law.

SORIN GROUP ITALIA warrants that all reasonable care has been taken in the manufacture of this medical device, as required by the nature of the device and the

use for which the device is intended.

SORIN GROUP ITALIA warrants that the medical device is capable of functioning as indicated in the current instructions for use when used in accordance with them by a qualified user and before any expiry date indicated on the packaging.

However, SORIN GROUP ITALIA cannot guarantee that the user will use the device correctly, nor that the incorrect diagnosis or therapy and/or that the particular physical and biological characteristics of an individual patient, do not affect the performance and effectiveness of the device with damaging consequences for the patient, even though the specified instructions for use have been respected.

SORIN GROUP ITALIA, whilst emphasizing the need to adhere strictly to the instructions for use and to adopt all the precautions necessary for the correct use of the device, cannot assume any responsibility for any loss, damage, expense, incidents or consequences arising directly or indirectly from the improper use of this device.

SORIN GROUP ITALIA undertakes to replace the medical device in the event that it is defective at the time of placing on the market or whilst being shipped by SORIN GROUP ITALIA up to the time of delivery to the final user unless such defect has been caused by mishandling by the purchaser.

The above replaces all other warranties explicit or implicit, written or verbal, including warranties of merchantability and fitness for purpose. No person, including any representative, agent, dealer, distributor or intermediary of SORIN GROUP ITALIA or any other industrial or commercial organization is authorized to make any representation or warranty concerning this medical device except as expressly stated herein. SORIN GROUP ITALIA disclaims any warranty of merchantability and any warranty of fitness for purpose with regard to this product other than what is expressly stated herein. The purchaser undertakes to comply with the terms of this Limited Warranty and in particular agrees, in the event of a dispute or litigation with SORIN GROUP ITALIA, not to make claims based on alleged or proven changes or alterations made to this Limited Warranty by any representative, agent, dealer, distributor or other intermediary.

The existing relations between the parties to the contract (also in the case that it is not drawn up in writing) to whom this Warranty is given as well as every dispute related to it or in any way connected to it as well as anything related to it or any dispute concerning this Warranty, its interpretation and execution, nothing excluded and/or reserved, are regulated exclusively by the Italian law and jurisdiction. The court chosen is the Court of Modena (Italy).

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