D 905 AVANTE

dideco

INSTRUCTIONS FOR USE

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MODE D'EMPLOI

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
INSTRUCTIONS FOR USE

I. CONTENTS

- Contents
- Description
- Technical features
- Intended use
- Safety information
- Set up
- Priming and recirculation procedure
- Initiating bypass
- Operation during bypass
- Terminating bypass
- Blood recovery after bypass
- Use of active venous drainage with vacuum
- Use of cardiotomy for post-operative autotransfusion
- Oxygenator replacement
- Medical devices for use with the D 903 AVANT
- Return of used products
- Limited Warranty

A. DESCRIPTION

The D 903 AVANT is a microporous hollow fibre membrane oxygenator consisting of a gas exchange module with an integrated heat exchanger and a Hardshell Reservoir divided into two compartments:

1. Venous Reservoir (fig. 2, ref.3)
2. Cardiotomy Reservoir with integrated cardiotomy filter (fig. 2, ref.4)

These two compartments are connected by means of a connection key (fig. 2, ref. 5). The device is single use, non-toxic, non pyrogenic, supplied **STERILE** and packaged singly. Sterilised by ethylene oxide. The level of ethylene oxide residuals in the device is within the limits established by national regulations in the country of use. The device is available in the following version:

- [A] D 903 AVANT with safety valve (30 µm cardiotomy filter)
[B] D 903 AVANT without safety valve (30 µm cardiotomy filter)

B. TECHNICAL FEATURES

- | | |
|---|---------------------------|
| - Recommended max blood flow | 7500 ml/min |
| - Reference blood flow (AAMI Standard) | 10000 ml/min |
| - Membrane type | Microporous Polypropylene |
| - Membrane surface area | 1.7 m ² |
| - Heat exchanger surface area | 0.14 m ² |
| - Hardshell Reservoir Capacity | 4700 ml |
| - Venous Reservoir | 3500 |
| - Cardiotomy Reservoir | 1200 |
| - Static priming volume (oxygenating module + heat exchanger) | 250 ml |
| - Connections: | |
| - Venous reservoir return | 1/2" (12.7 mm) |
| - Venous reservoir outlet | 3/8" (9.53 mm) |
| - Oxygenator venous inlet | 3/8" |
| - Oxygenating module arterial outlet | 3/8" |

C. INTENDED USE

The D 903 AVANT is intended for use in cardiopulmonary bypass circuits as a device to replace the function of the lungs, (transfer of oxygen and removal of carbon dioxide) to control the arterial/venous temperature and as a venous blood reservoir. The blood to be treated should contain anticoagulant. The D 903 AVANT is an ADULT oxygenator intended for use in procedures on patients with a body weight of less than 130 Kg (286 lb.). The D 903 AVANT should not be used longer than 6 hours. Contact with blood for longer periods is not advised. The D 903 AVANT should be used in combination with medical devices listed in section M (Medical devices for use with the D 903 AVANT).

D. SAFETY INFORMATION

Information intended to attract the attention of the user to potentially dangerous situations and to ensure correct and safe use of the device is indicated in the text in the following way:

WARNING

Indicates serious adverse reactions and potential safety hazards for practitioner and/or patient that can occur in the proper use or misuse of the device and also the limitations of use and the measures to be adopted in such cases.

CAUTION

Indicates any special care to be exercised by a practitioner for the safe and effective use of the device.

EXPLANATION OF THE SYMBOLS USED ON THE LABELS

For single use only (Do not reuse)

Batch code (number)
(reference for product traceability)

Use by (Expiry date)

Date of manufacture

Sterile - Ethylene oxide sterilised

Non Pyrogenic

Contains PHTHALATE

Latex free

Warning: Do not resterilize.

Contents sterile only if package is not opened, damaged or broken

Catalogue (code) number

Attention, see instruction for use

This way up

Fragile; handle with care

Quantity

Keep away from heat

Keep dry

The following is general safety information with the aim of advising the operator in preparing to use the device. Also, specific safety information is given in the instructions for use at locations in the text where that information is relevant for correct operation.

WARNING

- The device must be used in accordance with the instructions for use provided in this manual.
- The device is intended to be used by professionally trained personnel.
- SORIN GROUP ITALIA is not responsible for problems arising from inexperience or improper use.
- FRAGILE, handle with care.
- Do not expose to temperatures below 0°C (32°F) or above 60°C (140°F).
- Keep dry.
- Always give and maintain a correct dose and accurate monitoring

of the anticoagulant before, during and after the bypass.

- For single use and for single-patient use only. During use the device is in contact with human blood, body fluids, liquids or gases for the purpose of eventual infusion, administration or introduction into the body, and due to its specific design it cannot be fully cleaned and disinfected after use. Therefore, reuse on other patients might cause cross-contamination, infection and sepsis. In addition, the reuse increases the probability of product failure (integrity, functionality and clinical effectiveness).
- The device contains phthalates. Considering the nature of contact with the body, the limited duration of contact and the number of treatments per patient, the amount of phthalates which might be released from the device does not raise specific concerns about residual risks. Further information is available on request from Sorin Group Italia.
- The device must not undergo any further processing.
- Do not resterilise.
- After use, dispose of the device in accordance with applicable regulations in force in the country of use.
- The device must only be used if sterile. In the event of the device being supplied NOT STERILE (in which case, this is indicated by the wording "NOT STERILE" on the pack) contact SORIN GROUP ITALIA or an authorised dealer to agree on a sterilisation method.
- In the version with safety valve it is inadvisable to apply the vacuum to the cardiomy/venous reservoir.
- For further information and/or in case of complaint contact SORIN GROUP ITALIA or the authorised local representative.

CAUTION

Federal law (U.S.A.) restricts this device to sale by or on the order of a physician.

E. SET UP

1) POSITION THE HOLDER

Position the D 903 AVANT holder D 630 on the pump structure by means of the clamp at the upper end of the arm (fig.1, ref. 1).

2) FIX THE OXYGENATOR TO THE HOLDER

WARNING

- Sterility is guaranteed only if the sterile packaging is not wet, opened, damaged or broken. Do not use the device if sterility cannot be guaranteed.
 - Check the expiry date on the label attached. Do not use the device after the date shown.
 - The device must be used immediately after opening the sterile packaging.
 - The device must be handled aseptically.
- Remove the device from the sterile packaging.

WARNING

- Carry out a visual inspection and carefully check the device before use. Transport and/or storage conditions other than those prescribed may have caused damage to the device.
- Do not use solvents such as alcohol, ether, acetone, etc.: as contact may cause damage to the device.
- Do not allow halogenated liquids such as Halothane and Fluothane to come into contact with the polycarbonate housing of the device. This could cause damage which may compromise the integrity and proper functioning of the device.

Fix the D 903 AVANT onto the holder. The word "OPEN" on the water connectors locking system must be visible. Check that the notches on the water connectors locking system are aligned. Only then you can fix the oxygenator onto the holder. Insert the Hansen connectors and push D 903 AVANT down to the holder and turn the locking lever on "CLOSED" position. The D 903 AVANT will be correctly positioned only when the locking lever shows "CLOSED".

3) THERMOCIRCULATOR SET UP

Connect the water tubes to the holder by means of the female Hansen connectors SORIN GROUP ITALIA code 09028.

WARNING

- The use of different connectors from those indicated may cause resistance inside the circuit and reduce the efficiency of the heat exchanger.
- Do not obstruct the hole on the heat exchanger lower cover situated near the venous temperature probe. This hole connects the safety channel which helps prevent fluids crossing from one compartment to the other.
- The water temperature at the heat exchanger inlet must not exceed 42°C (108°F).
- The water pressure in the heat exchanger must not exceed 300 KPa (3 bar / 44 psi).

4) CHECK THE HEAT EXCHANGER

WARNING

Check the heat exchanger by recirculating water inside the heat exchanger for a few minutes. The integrity of the unit is guaranteed if there are absolutely no leaks from the water compartment or from the safety channel hole.

5) CIRCUIT CONNECTIONS

WARNING

All connections downstream of the pump must be secured by means of ties.

VENOUS LINE: Connect a venous line of 1/2" to the connector indicated on the venous reservoir as "VENOUS RETURN" (fig. 2, ref.15). The Venous Return connector can be rotated 180° to find the most convenient position of the venous tubing. This rotation is possible only after removing the polycarbonate insert that fixes the connector.

CAUTION

Remove the plastic insert only after connecting the venous tubing to the venous return connector.

ASPIRATION LINES:

Without separation of blood recovered by the aspirators:

after removing protective caps from the "filtered" inlet connectors on the top of the Cardiomy Reservoir (four 1/4" / 6.35 mm connectors, two 3/8" / 9.53 mm connectors), connect the aspiration lines and rotate the turrets (fig. 2, ref. 8) towards the suction pumps.

With separation of blood recovered by the aspirators:

connect the two "Y" adapters found in the convenience kit to the two connectors on the top of the Cardiomy ("Cardiomy Bypass Port"s (fig. 2, ref. 16)). Connect to the inlets of the two adapters (four 1/4" / 6.35 mm inlets, two Luer lock inlets) any lines that should connect directly to the Venous Reservoir: ventricular aspirators, aortic root aspirator, arterial filter drainage, blood concentrator, Venous Reservoir quick-priming. Finally, connect the remaining suction lines to the filtered inlets on the Cardiomy Reservoir (fig. 2, ref. 8). During the entire procedure, keep the connection (fig. 2, ref. 5) key closed (positioned at the bottom).

ARTERIAL LINE: remove the red cap on the oxygenator arterial outlet indicated as "ARTERIAL OUTLET" (fig. 2, ref. 7) and connect a 3/8" line.

PUMP LINE: the pump segment should be set up between the Venous Reservoir outlet connector (fig. 2, ref. 9) and the oxygenator venous inlet connector (fig. 2, ref. 6) taking account of the direction of rotation of the pump.

Remove yellow cap from the gas scavenging connector (fig. 2, ref. 18) of the Cardiomy Reservoir.

CAUTION

If oxygenated blood is necessary for blood cardioplegia, remove the red pos lock and connect the 1/4" blood line of the cardioplegia circuit to the D 903 AVANT coronary outlet port using the D523C reducer (provided with the product).

The coronary outlet port has a self-sealing valve which allows connection of the D 523C reducer also during extracorporeal

bypass without any leakage or spillage of fluid.

WARNING

If a connection is made to the coronary outlet port during bypass, the line to be connected must be unclamped and unpressurised such that blood will flow into it upon connection.

6) SAMPLING MANIFOLD

When sampling is needed from a more remote position remove the A/V manifold (fig. 2, ref. 10) from its housing and insert the plate in the "D 712 Holder for Sampling Manifold". The tube coil attached to the manifold allows it to be repositioned to a range of approximately 1 metre. Verify the connection of the manifold purge line with the "filtered" luer lock of the Venous Reservoir.

Check that the selector switch is in the "OFF" position.

7) CONNECT THE TEMPERATURE PROBES

The connection for the arterial temperature probe (red - fig. 2, ref. 11) is positioned next to the arterial outlet, while the venous probe site (blue - fig. 2, ref. 12) is at the side of the venous inlet.

SORIN GROUP ITALIA temperature probes have the code 09026.

8) CLOSE THE PURGING/RECIRCULATION LINE

Close the purging/recirculation stopcock (refer to the diagram on the label).

9) CONNECT THE GAS LINE

Remove the green cap from the gas inlet connector indicated as "GAS INLET" (fig. 2, ref. 13) and connect the 1/4" gas line. Ensure that the gas supply is from a suitable air/oxygen mixer such as the Sechrist, code 09046 (available from SORIN GROUP ITALIA) or a system with compatible technical features. A capnograph connector can be found in the centre of "GAS ESCAPE" (fig. 2, ref. 14) connector.

CAUTION

- The "GAS ESCAPE" system is designed to avoid any possible risk of blocking the gas outlet; such blockage could cause the immediate passage of air to the blood compartment.
- SORIN GROUP ITALIA recommends the use of a bubble trap or filter on the arterial line to reduce the risk of emboli transmission to the patient.

10) VAPOROUS ANAESTHETICS

The oxygenator is suitable for use with volatile anaesthetic isoflurane and sevoflurane, by mean of a suitable narcosis gas evaporator. If these vaporous anaesthetics are used, some method of scavenging the gas from the oxygenator should be considered.

The protocol, the concentration and the monitoring of the anaesthetic gases administrated to the patient, is under the sole responsibility of the physician in charge of the treatment.

WARNING

The only volatile anaesthetics suitable for this use are isoflurane and sevoflurane

WARNING

The methods adopted for vaporous anaesthetic gas scavenging should not increase or reduce in any way the pressure level at the oxygenator fibres.

F. PRIMING AND RECIRCULATION PROCEDURE

WARNING

Do not use alcoholic priming solutions: such solutions could compromise the proper function of the oxygenating module.

1) KEEP THE GAS FLOW OFF

2) KEEP THE OXYGENATOR PURGING/RECIRCULATION LINE CLOSED

Check that the purging/recirculation stopcock is closed.

3) CLOSE VENOUS AND ARTERIAL LINES

Clamp the venous line. Clamp the arterial line some centimetres away from arterial outlet connector of the oxygenator.

4) CHECK THE HEAT EXCHANGER

Verify again the integrity of the heat exchanger, paying particular attention to possible water leaks.

5) VENOUS RESERVOIR PRIMING

Secure with ties all aspiration lines connected to the Cardiotomy reservoir. Fill the Cardiotomy Reservoir with sufficient liquid to ensure the intended haematocrit is obtained, taking into account:

- the static priming volume of the oxygenator is 250 ml;
- the 3/8" tube capacity is 72 ml/m;
- the 1/2" tube capacity is 127 ml/m.

Clamp the Venous Reservoir outlet.

In order to fill the Venous Reservoir or if the Cardiotomy Reservoir capacity (1200 mls) is not enough, open the connection to the Venous Reservoir by raising the connection key (fig. 2, ref.5) on the top of the Cardiotomy Reservoir.

6) CIRCUIT PRIMING

WARNING

- The pressure level inside the blood compartment of the oxygenating module shall not exceed 100 KPa (1 bar / 14 psi).

Clamp the venous reservoir outlet and remove the pump segment from the arterial pump head. Fill the pump segment by keeping it at the same height as the Venous Reservoir and slowly opening the clamp occluding it. By slowly bending downward the tube segment to be filled, the air in the tube will be routed to the oxygenator. The complete priming of the oxygenating module is completed by gravity. Once the D 903 AVANT is filled, place the pump segment in the arterial pump.

7) OPEN VENOUS AND ARTERIAL LINES

Remove the clamp from the venous and arterial lines and increase flow up to 2000 ml/min.

8) OPEN THE PURGING/RECIRCULATION LINE

Once the steps up to and including point 7 have been carried out, increase the arterial pump speed until the flow reaches the maximum value of 7.5 l/min. Open the purging/recirculation stopcock for some seconds in order to prime the oxygenator purging/recirculation line.

9) PURGE THE AIR CONTAINED IN THE CIRCUIT

During this phase it is necessary to tap the entire circuit in order to facilitate the removal of microbubbles from the tube walls. After some minutes in which the flow is maintained at a high rate, all air will be evacuated.

10) PRIME OF THE SAMPLING MANIFOLD

The removal of the air from the A/V sampling system and from the coronary outlet port is achieved as follows:

- Rotate the manifold selector switch in the "A-V SHUNT" position.
- Check that the A/V sampling line is completely purged.
- Reset the manifold selector switch in the OFF position.

11) CLOSE THE PURGING/RECIRCULATION LINE

After 3-5 minutes in which the flow is maintained at a high rate, all air will be evacuated and it is possible to close the purging/recirculation line from the dedicated stopcock.

12) CLOSE THE VENOUS AND ARTERIAL LINES

CAUTION

- Do not use pulsatile flow during priming.
- Check the correct dosage of anticoagulant in the system before starting the bypass.
- SORIN GROUP ITALIA recommends the use of the pump speed control to reduce or stop the arterial flow slowly.
- Do not use the pump on/off switch until the pump speed is zero.

WARNING

- If the reduction connector (D523C) and a circuit have been connected to the coronary outlet port, check the priming of the connected line.
- Clamp the line some centimetres away from the outlet.
- Do not create a negative pressure at the coronary outlet. Negative pressure in the blood compartment could cause microbubble formation.

G. INITIATING BYPASS

1) OPEN THE ARTERIAL AND VENOUS LINES

Remove first the clamp from the arterial line, then remove the clamp on the venous line. Start bypass with a blood flow appropriate to patient size.

Check constantly the blood level in the Venous Reservoir.

2) CHECK THE CORRECT OPERATION OF THE HEAT EXCHANGER

Check the temperature of the venous and arterial blood.

3) SELECTION OF THE APPROPRIATE GAS FLOW

The suggested gas/blood flow ratio in normothermia is 1:1 with a FiO_2 of 80:100%.

WARNING

- Always open the gas flow after the blood flow. The gas/blood flow ratio must never exceed 2:1.
- The pressure in the blood compartment must always exceed that of the gas compartment. This is to prevent gas emboli appearing in the blood compartment.

4) BLOOD GAS MONITORING

After a few minutes of bypass operation, measure the gas content of the blood. Depending on the values found, adjust the relevant parameters as follows:

High pO_2	→	Decrease FiO_2
Low pO_2	→	Increase FiO_2
High pCO_2	→	Increase gas flow
Low pCO_2	→	Decrease gas flow

H. DURING BYPASS

1) CHECK THE VENOUS RETURN

If a higher venous return flow is necessary lower both the oxygenator and the venous reservoir with respect to the patient position.

WARNING

- The ACT (Activated Coagulation Time) must always be longer than or equal to 480 seconds in order to ensure adequate anticoagulation of extracorporeal circuit.
- If administration of anticoagulant to the patient is envisaged, use the "DRUGS PORT" of the sampling stopcock.

2) ARTERIAL SAMPLING

- Turn the selector switch on the sampling manifold (fig. 2, ref. 10) to the "A-SAMPLE PURGE" position (the arterial line will then be automatically purged).
- Insert the syringe into the sampling female luer marked SAMPLE PORT.
- Aspirate at least 2 ml of blood (the automatic purge avoids the need to repeat the operation). A self-closing valve will avoid any blood leakage from the luer.

3) VENOUS SAMPLING

- Turn the selector switch on the sampling stopcock to the "V-SAMPLE PURGE" position.
- Insert the syringe into the female luer marked "SAMPLE PORT".
- Aspirate and purge the blood sample a couple of times before proceeding to the venous sampling analysis. This double operation (aspiration and purge) can be done whilst maintaining the selector switch in the same position.

4) DRUGS INJECTION

- Turn the selector switch on the sampling stopcock to the "DRUGS INJECTION" position.
- Insert the syringe containing the drug into the female luer marked "DRUGS PORT" and inject the liquid.
- Turn the switch to the A-V SHUNT position. The system will be automatically flushed and the drug will flow into the venous line.

5) LOW FLOW RECIRCULATION

(Hypothermia associated with circulatory arrest).

- Reduce the gas flow to less than 500 ml/min.
- Open the recirculation line (lever of the purging/recirculation stopcock in the "RECIRCULATION" position) and clamp the Venous Reservoir inlet line (fig. 2, ref. 15).
- Reduce the flow from the arterial pump to 2000 ml/min.
- Clamp the oxygenator arterial line (fig. 2, ref. 7).

- Recirculate at a maximum flow of 2000 ml/min. throughout the patient's circulatory arrest.
- To restart bypass after circulatory arrest, open the venous and arterial lines and slowly increase the blood flow.
- Close the recirculation line (lever in "CLOSE" position).
- Adjust gas flow.

6) VENOUS/CARDIOTOMY RESERVOIR MANAGEMENT

As reported in the Description (chapter A), the D903 AVANT enables the separation of the Cardiomy Reservoir from the Venous Reservoir. This separation is achieved by opening (up position) or closing (down position) the connection lever located on the top of the Cardiomy Reservoir.

This option is suggested when the liquid from the suckers must be separated from the venous blood in order to be:

- rejected
- processed apart

If you need to add priming solution directly into the Venous Reservoir use the indicated port "Cardiotomy Bypass Port" (fig. 2, ref. 16). With the connection key in up position the hardshell reservoir acts as a standard Venous Reservoir.

WARNING

FOR DEVICES EQUIPPED WITH SELF-CLOSING SAFETY VALVE:

- The self-closing safety valve is a safety device to prevent massive air transmission to the patient in case the venous reservoir suddenly empties. The valve works with flows between 1 to 7.5 l/min. However, the self-closing safety valve does not relieve the perfusionist from carefully checking the level in the Venous Reservoir.

More attention has to be paid when emptying the reservoir with low blood level and/or low flows.

- In the case of insufficient venous return flow the self-closing safety valve might collapse. If this happens, consequent negative pressure might cause gas extraction from the blood if the arterial pump is not stopped immediately. To eliminate air from the oxygenating module, proceed as follows:

- Turn the gas flow off.
- Turn the arterial pump off.
- Clamp the arterial line.
- Restore the required volume of liquid into the Venous Reservoir.
- Clamp the venous line.
- Open the purging/recirculation line and recirculate at a rate of 2000 ml/min until air removal from the system is completed.
- Restart bypass by opening the arterial and venous lines.
- Close the purging/recirculation line.
- Open gas flow.

7) CONTINUOUS AIR PURGE

The purging/recirculation stopcock in "PURGE" position features the continuous air purge diverting from the gas exchange module during bypass. In this condition and at full arterial blood flow, the continuous purge diverts from the arterial line only few mls/min.

I. TERMINATING BYPASS

Must be carried out after consideration of each individual patient's state. Act as follows:

- Turn the gas flow off.
- Turn the thermocirculator off.
- Slowly decrease the arterial flow to zero while closing the venous line.
- Clamp the arterial line.
- Open the recirculation line.
- Increase arterial flow until 2000 ml/min.
- In the case of separation of blood recovered by the aspirators:
 - remove the pos-lock on the connection key (fig. 2, ref. 5) and connect the adapter D 523C (provided with the product)
 - recover the blood collected in the Cardiomy Reservoir with an autotransfusion machine connected to the adapter with a 1/4" line
- wash the collected blood and reinfuse to the patient

CAUTION

- If extracorporeal circulation has to be subsequently restarted, a

minimum blood flow inside the D 903 AVANT must be maintained (maximum 2000 ml/min).

- During recirculation do not turn the thermocirculator off.
- Verify that the cardioplegia circuit connected to the coronary outlet port is properly clamped.

J. BLOOD RECOVERY AFTER BYPASS

- 1) Recover into the Venous Reservoir as much blood as possible from the venous line, as soon as the surgeon has removed the cannulae from the patient's vena cava.
- 2) Deliver blood into the aortic canulae as required by the patient's condition, slowly decreasing the level in the Venous Reservoir.
- 3) When the reservoir is nearly empty stop the arterial pump and clamp the arterial line.

K. USE OF ACTIVE VENOUS DRAINAGE WITH VACUUM

This method may be applied at any time during extra-corporeal circulation, provided that the steps below are followed. Using the kit code 086000, or equivalent supplied separately, and a vacuum regulation device, AVANT can be used for active venous drainage using vacuum. This technique constitutes an alternative to venous drainage by gravity and allows the use of shorter venous lines with reduced diameter and smaller-gauged cannulas.

L. USE OF CARDIOTOMY RESERVOIR FOR POST-OPERATIVE AUTOTRANSFUSION

If the use of the Cardiomy Reservoir for post-operative autotransfusion is envisaged, act as follows:

1. Disconnect the purging/recirculation line.
2. Separate the venous reservoir from the oxygenating module by removing the white hook.
3. Locate the venous reservoir on the "Post Operative Chest DrainageHolder code 05039" and use one of the following optional kits referring to their respective instructions for use:
 - D 540 AUTOTRANSFUSION CONVERSION KIT code 05053;
 - D 540 W AUTOTRANSFUSION CONVERSION KIT with water seal, code 05062.
4. If wishing to use the full capacity of the reservoir, lock the connection key (fig. 2, ref. 5) with the red clamp found in the convenience kit provided with the product.

M. OXYGENATOR REPLACEMENT

A spare oxygenator must always be available during perfusion. After 6 hours of use with blood or if particular situations occur, which may lead the person responsible for perfusion to determine that the safety of the patient may be compromised, (insufficient oxygenator performance, leaks, abnormal blood parameters etc.), proceed as follows for oxygenator replacement.

CAUTION

Use sterile methods during all replacement procedure.

Replacement of the oxygenator and the AVANT DUALRESERVOIR

- 1) Turn the gas flow off.
- 2) Close, by means of a double clamp, the venous line (5 centimetres apart).
- 3) Turn the arterial pump off and close, by means of a double clamp, the arterial line (5 centimetres apart) placed next to the oxygenator.
- 4) Turn the thermocirculator off, clamp and remove the water lines.
- 5) Disconnect the gas line, all monitoring and sampling lines.
- 6) Cut the venous return and the arterial lines in the section between the two clamps, leaving a sufficient length of tubing to allow re-connection.
- 7) Remove the D 903 AVANT from the holder and remove the pump segment from the arterial pump.
- 8) Place a new D 903 AVANT on the holder. Connect all lines (i.e. venous to the Venous Reservoir, arterial and gas to the oxygenator, pump line to Venous Reservoir and oxygenator).

WARNING

In this phase, keep the venous and arterial lines clamped.

- 9) Open the water lines on the holder, turn the thermocirculator on

and check the integrity of the new D 903 AVANT.

- 10) Fill the cardiomy reservoir of the new D 903 AVANT with priming solution through the 1/4" or 3/8" quick prime connectors located on the head of the venous reservoir.
- 11) Prime the new D 903 AVANT and evacuate the microbubbles, as described in the priming and recirculation procedure.
- 12) Verify all connections and secure with ties.
- 13) Remove clamps from the venous and the arterial line, close the purging/recirculation line and start the bypass again.
- 14) The blood remaining in the replaced Venous Reservoir may be recovered by connecting its outlet port to one of the 3/8" inlet connectors of the new reservoir.
- 15) The blood contained in the oxygenator and heat exchanger may be poured into the new Venous Reservoir by connecting the arterial line to one of the 3/8" inlet connectors of the new reservoir.

Replacement of only the oxygenating module

- 1) Turn the gas flow off and disconnect the gas line.
- 2) Clamp the venous return.
- 3) Turn the arterial pump off and close, by means of a double clamp, the arterial line (5 centimetres apart) next to the oxygenator.
- 4) Close, by means of a double clamp, the oxygenator inlet line near the oxygenating module connector (5 centimetres apart).
- 5) Turn the thermocirculator off, clamp and remove the water lines.
- 6) Close the purging/recirculation stopcock and disconnect the purging/recirculation line.
- 7) Clamp and disconnect the cardioplegia line, if connected.
- 8) Turn the sampling manifold selector switch to the "OFF" position.
- 9) Disconnect the arterial sampling avoiding any contamination of the sampling stopcock luer lock (i.e. by connecting it to a female luer lock placed on the top of the reservoir).
- 10) Remove all other monitoring and sampling lines from the oxygenating module.
- 11) Cut the oxygenator inlet line and the arterial line in the section between the two clamps, leaving a sufficient length of tubing to allow re-connection.
- 12) Move the water locking system to the "OPEN" position.
- 13) Lift the oxygenator and rotate the hook (fig. 1, ref. 2) on the upper part of the holder.
- 14) Fit the slot of the reservoir into the hook of the holder.
- 15) Uncouple the reservoir from the oxygenating module by means of the white interface (fig. 2, ref.19)
- 16) Remove the oxygenating module from the holder.
- 17) Place a new D 903 AVANT oxygenating module on the holder.
- 18) Fix it by moving the water connector locking system to the "CLOSE" position.
- 19) Connect the oxygenator inlet line and the arterial line to the oxygenating module

In this phase, keep the venous and arterial lines clamped.

- 20) Switch on the thermocirculator and check the heat exchanger.
- 21) Lift the reservoir and rotate the hook. Put the reservoir on the oxygenating module. If possible couple the reservoir and the oxygenating module by means of the white interface.
- 22) Connect the gas line, the arterial sampling line and all the monitoring lines.
- 23) Remove the clamp from the venous line.
- 24) Keeping the recirculation line open (purging/recirculation stopcock in "RECIRCULATION" position) prime the oxygenating module through the arterial pump. Purge the air contained in the module at a maximum flow rate of 2000 ml/min
- 25) Close the recirculation line (purging/recirculation stopcock in "CLOSE" position).
- 26) Remove the clamp from the arterial line and restart the bypass.
- 27) Connect, if necessary, the cardioplegia line and purge it.

N. MEDICAL DEVICES FOR USE WITH THE D 903 AVANT

For post-operative autotransfusion with the Venous Reservoir one of the two following kits should be used:

- D 540 AUTOTRANSFUSION CONVERSION KIT code 05053.
- D 540 W AUTOTRANSFUSION CONVERSION KIT with water seal, code 05062.

All tubing used to make the circuit connections must be of a diameter

which is compatible with the dimensions of the connectors on the device (3/8", 1/4", 1/2"). Temperature controls must be carried out using SORIN GROUP ITALIA probes, code 09026 compatible to YSI Series 400. Use Sechrist air/oxygen mixer (SORIN GROUP ITALIA code 09046) or a system with compatible technical features. Any heating/cooling system (thermocirculator) may be used, provided that the connectors to the water distributor holder are of the Hansen type (SORIN GROUP ITALIA code 09028). Currently SORIN GROUP ITALIA is not aware of any contraindications to the use of the device with occlusive or non-occlusive peristaltic pumps or with centrifugal pumps. The use of other types of pump must be agreed with SORIN GROUP ITALIA.

O. RETURN OF USED PRODUCTS

Should the user be dissatisfied with anything related to the quality of the product, the product distributor or the authorized local SORIN GROUP ITALIA representative should be notified.

All parameters considered critical by the user must be reported with particular care and urgency. The following is the minimum information that should be provided:

- Detailed description of the event and, if pertinent, the conditions of the patient;
- Identification of the product involved;
- Lot number of the product involved;
- Availability of the product involved;
- All the indications the user considers useful in order to understand the origin of the elements of dissatisfaction.

SORIN GROUP ITALIA reserves the right to authorize, if necessary, recall of the product involved in the notification for assessment. If the product to be returned is contaminated, it must be treated, packed and handled in conformity with the provisions of the legislation in force in the country where the product was used.

It is the responsibility of the health care institution to adequately prepare and identify the product for return shipment. Do not return products that have been exposed to blood borne infectious diseases.

ONLY for US customers

If for any reason the product must be returned to the manufacturer, a returned good authorisation (RGA) number is required from Sorin Group USA cardiovascular Inc. prior to shipping.

If the product has been in contact with blood or blood fluids, it must be thoroughly cleaned and disinfected before packing. It should be shipped in either the original carton or an equivalent carton to prevent damage during shipment, and it should be properly labelled with an RGA number and an indication of the biohazardous nature of the content in the shipment.

The shipping address for returned goods in the US is:

Sorin Group USA, Inc.
Returned CV Products
14401 West 65th Way
Arvada, CO 80004-3599
FAX (800) 323 4031.

P. LIMITED WARRANTY

This Limited Warranty is in addition to any statutory rights of the Purchaser pursuant to applicable law.

SORIN GROUP ITALIA warrants that all reasonable care has been taken in the manufacture of this medical device, as required by the nature of the device and the use for which the device is intended.

SORIN GROUP ITALIA warrants that the medical device is capable of functioning as indicated in the current instructions for use when used in accordance with them by a qualified user and before any expiry date indicated on the packaging. However, SORIN GROUP ITALIA cannot guarantee that the user will use the device correctly, nor that the incorrect diagnosis or therapy and/or that the particular physical and biological characteristics of an individual patient, do not affect the performance and effectiveness of the device with damaging conse-

quences for the patient, even though the specified instructions for use have been respected. SORIN GROUP ITALIA, whilst emphasizing the need to adhere strictly to the instructions for use and to adopt all the precautions necessary for the correct use of the device, cannot assume any responsibility for any loss, damage, expense, incidents or consequences arising directly or indirectly from the improper use of this device. SORIN GROUP ITALIA undertakes to replace the medical device in the event that it is defective at the time of placing on the market or whilst being shipped by SORIN GROUP ITALIA up to the time of delivery to the final user unless such defect has been caused by mishandling by the purchaser. The above replaces all other warranties explicit or implicit, written or verbal, including warranties of merchantability and fitness for purpose. No person, including any representative, agent, dealer, distributor or intermediary of SORIN GROUP ITALIA or any other industrial or commercial organization is authorized to make any representation or warranty concerning this medical device except as expressly stated herein. SORIN GROUP ITALIA disclaims any warranty of merchantability and any warranty of fitness for purpose with regard to this product other than what is expressly stated herein. The purchaser undertakes to comply with the terms of this Limited Warranty and in particular agrees, in the event of a dispute or litigation with SORIN GROUP ITALIA, not to make claims based on alleged or proven changes or alterations made to this Limited Warranty by any representative, agent, dealer, distributor or other intermediary. The existing relations between the parties to the contract (also in the case that it is not drawn up in writing) to whom this Warranty is given as well as every dispute related to it or in any way connected to it as well as anything related to it or any dispute concerning this Warranty, its interpretation and execution, nothing excluded and/or reserved, are regulated exclusively by the Italian law and jurisdiction.

The court chosen is the Court of Modena (Italy).

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